

**IN THE UNITED STATES PATENT AND TRADEMARK OFFICE**

Application of: DeStefano, G. et al. ) Confirm. No.: 7531  
Serial No.: 10/804,710 ) Art Unit: 1616  
Filed: March 19, 2004 ) Examiner: Haghighatian, Mina  
For: Formulation for a metered dose inhaler using hydro-fluoro-alkanes  
as propellants  
Docket No.: 9/277 (539/168)

Mail Stop Amendment  
Commissioner for Patents  
P.O. Box 1450  
Alexandria, VA 22313-1450

DECLARATION OF GEORGE DESTEFANO UNDER 37 C.F.R. § 1.132

I, George DeStefano, declare:

1. This declaration is submitted to present evidence in response to the Office Action dated October 4, 2007 in the above-referenced patent application.

2. I am employed at Boehringer Ingelheim Pharmaceuticals, Inc. I am an inventor in the above-referenced patent application.


3. It is noted that the Office Action contains a rejection of claims 5 and 7 under 35 U.S.C. § 103(a) as being unpatentable over Lewis et al. (EP 1219293) in view of Jager et al. (WO 9413262), and of claims 1-7 on the ground of nonstatutory obviousness-type double patenting as being unpatentable over claims 1-11 of U.S. Pat. No. 6,423,298 in view of Lewis et al. The documents submitted herein in the attached pages provide additional experimental data showing the criticality of the claimed water content in the claimed formulations.

4. The following Analytical Reports provide support for the criticality of the water content in suspension formulations (Analytical Report, *Results for In Use 6 Hour Delay Single Actuation Reproducibility Study for Ipratropium Bromide/Albuterol Sulfate (HFA-134a) Inhalation Aerosol, with Water Contents Ranging from Inherent Water to 2500ppm*,

File No. AR-030012, February 5, 2003, Boehringer Ingelheim Pharmaceuticals Inc.; Analytical Report, *Evaluation of the Effect of Water on Single Actuation Reproducibility for Ipratropium Bromide/Albuterol Sulfate (HFA-134a) Inhalation Aerosol 0.021/0.120 mg TTV, 10 mL Equipped with the Bespak Valve Actuated using at Patient In-Use Testing Scheme*, Report No. AS/C-03004, January 28, 2003, Boehringer Ingelheim Pharmaceuticals Inc.). The data presented in Report No. AS/C-03004 show that a water concentration of 1500 ppm and greater (equivalent to 0.15% water content) in a suspension formulation containing albuterol sulphate and ipratropium bromide is needed to ensure that the single actuation reproducibility difference does not impact the product (see Abstract, Report No. AS/C-03004). File No. AR-030012 shows that water content of 300 to 1200 ppm (equivalent to 0.03 to 0.12% water content) in a suspension formulation containing albuterol sulphate and ipratropium bromide results in poor reproducibility, while a water content of 1500 to 2500 ppm (equivalent to 0.15 to 0.25% water content) exhibited good reproducibility (page 29, File No. AR-030012). This shows that the water content of the suspension formulation is a critical parameter in the single actuation reproducibility over time.

5. The Analytical Reports demonstrate the criticality of the water content in suspension formulations to achieve single actuation reproducibility of inhalants. The critical range of water content in a suspension formulation containing albuterol sulphate and ipratropium bromide to achieve reproducibility is from 0.13% and 0.25%.

6. I declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and that these statements were made with the knowledge that willful false statements and the likes so made are punishable by fine or imprisonment, or both under Section 1001 of Title 18 of the United States Code, and that such willful false statements may jeopardize the validity of the application or any patents issuing thereon.

  
George DeStefano  
Date: March 10, 2008



**Boehringer  
Ingelheim**

**Report Number:**

**AS/C-03004**

**ANALYTICAL SCIENCES DEPARTMENT  
CONTRACT LABORATORY REPORT**

**Title:** Evaluation of the Effect of Water on Single Actuation Reproducibility for Ipratropium Bromide/Albuterol Sulfate ( ) Inhalation Aerosol 0.021/0.120 mg TTV, 10 mL Equipped with the Bespak Valve Actuated using at Patient In-Use Testing Scheme

**Laboratory:** Cardinal Health RTP

**Report Date:** Jan. 28, 2003

**Product/  
Project:** Ipratropium Bromide/Albuterol Sulfate  
Inhalation Aerosol 0.021/0.120 mg TTV, 10 mL

**BI Project Code:**

**Person Submitting  
Report for Archival:** Paul D. Curry, Jr., Ph.D.

**Contract Lab  
Protocol #:**

**Date Report  
Archived:** 1/29/03

**Key Words:** Ipratropium Bromide/Albuterol Sulfate ( ) Inhalation Aerosol 0.021/0.120 mg TTV, 10 mL, Effect of Water, Single Actuation Reproducibility, Bespak Valve, Patient In-Use

**Abstract:** It is well known that additional water can have a detrimental effect on the performance of non-aqueous suspensions. This study was designed to study the impact of water on single actuation reproducibility (SAR) and to initially identify a minimum water concentration where the reproducibility is suitable. Other studies have shown low dose delivery for the first shot followed by a higher dose for the second for albuterol sulfate at low water concentrations. The addition of water to the formulation brings the dose for the individual actuation back to the target delivery. Ipratropium bromide does not show this behavior because it is in solution in this formulation unlike albuterol sulfate which is suspended.

This study demonstrates the same behavior. The data would indicate that water concentration of 1500 ppm and greater are needed to ensure that the single actuation reproducibility difference does not impact the product.

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7501

Theoretical Water Content: Inherent

Analyst: Brook White

| Actuation Number | Test | Inverted            |                    |         |       |         |       |         |       |         |       |
|------------------|------|---------------------|--------------------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1               |                    | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
| 25               | SAR  | Alb SO4             | IPB                | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 26               | SAR  | 115.13              | 20.56              | 118.21  | 20.85 | 115.17  | 20.64 | 117.57  | 20.63 | 122.36  | 21.72 |
| 33               | SAR  | 113.42              | 19.74              | 112.10  | 19.45 | 117.79  | 20.54 | 117.26  | 20.34 | 119.82  | 21.42 |
| 34               | SAR  | 77.03               | 19.66              | 66.59   | 19.61 | 90.17   | 19.95 | 62.84   | 19.43 | 82.38   | 19.42 |
| 45-46            | Dose | 245.57              | 20.14              | 202.70  | 19.47 | 216.22  | 19.59 | 238.14  | 21.24 | 214.94  | 20.25 |
| 69               | SAR  | 159.74              | 19.46              | 139.27  | 18.83 | 144.63  | 19.52 | 169.23  | 19.22 | 149.47  | 19.76 |
| 70               | SAR  | 60.95               | 19.29              | 75.53   | 19.87 | 76.67   | 19.34 | 70.93   | 18.94 | 128.19  | 18.60 |
| 105-106          | Dose | 240.53              | 20.28              | 237.87  | 19.34 | 256.88  | 20.88 | 232.40  | 20.84 | 222.82  | 21.00 |
| 129              | SAR  | 133.76              | 19.27              | 95.46   | 19.54 | 95.94   | 20.03 | 92.92   | 19.78 | 95.99   | 19.87 |
| 130              | SAR  | 169.42 <sup>1</sup> | 19.62 <sup>1</sup> | 57.22   | 19.31 | 79.18   | 19.93 | 69.29   | 18.84 | 56.62   | 18.79 |
| 165-166          | Dose | 161.63 <sup>2</sup> | 19.62 <sup>2</sup> | 181.90  | 19.69 | 156.13  | 20.43 | 177.50  | 20.78 | 161.25  | 19.90 |
| 187              | SAR  | 103.76              | 19.55              | 101.53  | 19.61 | 104.51  | 20.45 | 100.96  | 19.96 | 92.65   | 19.65 |
| 188              | SAR  | 87.58               | 19.52              | 58.20   | 25.05 | 47.18   | 19.52 | 57.62   | 19.52 | 92.57   | 19.22 |
|                  |      | 146.99              | 20.16              | 168.57  | 20.82 | 142.08  | 21.13 | 153.97  | 20.93 | 136.70  | 20.66 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 263   | 257   | 260   | 262   | 265   |
| 199-203   | 364   | 343   | 343   | 331   | 325   |

<sup>1</sup>These values were confirmed by reinjection since they were not in trend with the data set. The reinjected samples yielded results of 167.16 and 18.81 µg for Albuterol Sulfate and Ipratropium Bromide respectively.

<sup>2</sup>These values were confirmed by reinjection since they were not in trend with the data set. The reinjected samples yielded results of 160.32 and 19.42 µg for Albuterol Sulfate and Ipratropium Bromide respectively.

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7501

Theoretical Water Content: Inherent

Analyst: Nichole Bynum

| Actuation Number | Test | Upright |       |                     |                    |         |       |                    |                    |                    |                    |
|------------------|------|---------|-------|---------------------|--------------------|---------|-------|--------------------|--------------------|--------------------|--------------------|
|                  |      | Can 6   |       | Can 7               |                    | Can 8   |       | Can 9              |                    | Can 10             |                    |
|                  |      | Alb SO4 | IPB   | Alb SO4             | IPB                | Alb SO4 | IPB   | Alb SO4            | IPB                | Alb SO4            | IPB                |
| 25               | SAR  | 114.77  | 20.55 | 121.29              | 21.97              | 121.85  | 21.73 | 115.84             | 20.69              | 121.53             | 21.52              |
| 26               | SAR  | 119.57  | 20.87 | 118.26              | 20.88              | 117.92  | 22.26 | 121.20             | 21.46              | 119.29             | 21.36              |
| 33               | SAR  | 39.60   | 19.58 | 47.68               | 20.48              | 42.23   | 20.70 | 40.38              | 20.39              | 47.37              | 20.29              |
| 34               | SAR  | 164.26  | 19.83 | 155.10              | 19.75              | 156.82  | 19.18 | 169.58             | 19.66              | 156.61             | 19.84              |
| 45-46            | Dose | 98.77   | 19.56 | 99.56               | 19.85              | 99.49   | 20.45 | 104.71             | 20.18              | 98.57              | 20.05              |
| 69               | SAR  | 46.50   | 19.64 | 50.50               | 19.16              | 41.45   | 19.82 | 50.31              | 20.04              | 52.55              | 20.33              |
| 70               | SAR  | 155.93  | 19.89 | 157.65              | 20.68              | 159.42  | 19.32 | 150.98             | 18.73              | 158.94             | 20.84              |
| 105-106          | Dose | 88.03   | 18.82 | 96.15               | 20.00              | 94.30   | 19.59 | 92.94              | 19.42              | 106.61             | 20.44              |
| 129              | SAR  | 51.96   | 20.29 | 55.25               | 20.24              | 50.28   | 12.78 | 89.13              | 61.90              | 78.36              | 17.70              |
| 130              | SAR  | 150.76  | 20.17 | 141.95              | 19.57              | 142.01  | 19.74 | 130.40             | 18.61              | 136.65             | 20.80              |
| 165-166          | Dose | 94.23   | 19.57 | 102.94 <sup>2</sup> | 17.66 <sup>2</sup> | 101.73  | 18.17 | 99.60 <sup>2</sup> | 16.83 <sup>2</sup> | 98.36 <sup>2</sup> | 16.97 <sup>2</sup> |
| 187              | SAR  | 48.15   | 20.44 | 60.41               | 20.65              | 42.98   | 20.22 | 46.04              | 20.68              | 60.89              | 22.05              |
| 188              | SAR  | 149.23  | 20.66 | 140.19              | 20.92              | 150.85  | 21.02 | 152.92             | 20.66              | 138.40             | 21.78              |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 6 | Can 7 | Can 8 | Can 9 | Can 10           |
|-----------|-------|-------|-------|-------|------------------|
| 16-20     | 275   | 270   | 260   | 268   | 279              |
| 199-203   | 341   | 337   | 327   | 330   | 402 <sup>1</sup> |

<sup>1</sup>The analyst observed an unusually weak plume for all 5 actuations.

<sup>2</sup>The analyst observed that the first actuation sounded weaker than normal.

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7502

Theoretical Water Content: 800ppm

Analyst: Dan Milka

| Actuation Number | Test | Inverted |       |         |       |         |       |         |       |         |       |
|------------------|------|----------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1    |       | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
| 25               | SAR  | Alb SO4  | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 26               | SAR  | 116.85   | 20.57 | 118.69  | 20.61 | 116.57  | 20.81 | 117.72  | 20.74 | 112.33  | 19.75 |
| 33               | SAR  | 116.44   | 19.92 | 118.78  | 20.30 | 113.46  | 19.87 | 118.76  | 20.62 | 110.60  | 19.21 |
| 34               | SAR  | 79.59    | 20.67 | 92.27   | 20.84 | 66.80   | 20.18 | 80.28   | 20.41 | 69.10   | 20.04 |
| 45-46            | Dose | 145.92   | 19.89 | 140.63  | 19.98 | 151.68  | 20.06 | 143.55  | 20.54 | 145.53  | 19.59 |
| 69               | SAR  | 110.49   | 19.82 | 108.81  | 20.23 | 111.42  | 20.15 | 108.24  | 20.16 | 106.46  | 19.49 |
| 70               | SAR  | 78.14    | 20.17 | 70.83   | 20.56 | 108.78  | 18.75 | 96.90   | 19.75 | 64.81   | 19.97 |
| 105-106          | Dose | 160.24   | 20.01 | 152.88  | 20.09 | 140.17  | 19.52 | 151.44  | 19.79 | 147.23  | 18.80 |
| 129              | SAR  | 108.82   | 19.64 | 112.10  | 19.69 | 109.69  | 19.58 | 108.86  | 19.80 | 106.69  | 19.31 |
| 130              | SAR  | 68.67    | 20.47 | 73.18   | 20.45 | 77.48   | 20.10 | 72.12   | 20.29 | 99.65   | 18.57 |
| 165-166          | Dose | 156.79   | 19.58 | 167.49  | 20.08 | 148.12  | 19.88 | 154.79  | 19.99 | 134.07  | 18.99 |
| 187              | SAR  | 112.93   | 19.43 | 111.55  | 18.86 | 100.26  | 19.74 | 104.20  | 20.03 | 100.62  | 19.23 |
| 188              | SAR  | 67.20    | 21.10 | 103.12  | 19.25 | 95.01   | 20.05 | 95.15   | 19.83 | 92.73   | 19.36 |
|                  |      | 147.90   | 20.39 | 124.93  | 20.32 | 123.27  | 19.63 | 131.34  | 20.49 | 116.31  | 19.92 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 688   | 722   | 714   | 735   | 686   |
| 199-203   | 762   | 811   | 802   | 821   | 819   |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7502

Theoretical Water Content: 800ppm

Analyst: Kathryn Johnson

| Actuation Number | Test | Upright |       |         |       |         |       |         |       |         |       |
|------------------|------|---------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 6   |       | Can 7   |       | Can 8   |       | Can 9   |       | Can 10  |       |
|                  |      | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 115.28  | 20.70 | 119.93  | 21.20 | 115.46  | 20.36 | 113.94  | 20.15 | 119.07  | 21.20 |
| 26               | SAR  | 114.34  | 20.09 | 124.94  | 20.97 | 117.20  | 20.01 | 118.16  | 20.34 | 118.37  | 20.66 |
| 33               | SAR  | 39.99   | 19.47 | 41.84   | 20.47 | 46.49   | 19.44 | 38.20   | 21.41 | 43.36   | 20.82 |
| 34               | SAR  | 157.71  | 18.38 | 164.30  | 19.43 | 157.55  | 19.46 | 168.46  | 20.03 | 166.32  | 19.99 |
| 45-46            | Dose | 96.73   | 19.35 | 92.46   | 18.89 | 93.65   | 19.13 | 94.80   | 18.68 | 100.58  | 19.88 |
| 69               | SAR  | 38.59   | 19.78 | 35.71   | 19.36 | 49.57   | 19.41 | 38.83   | 19.89 | 42.23   | 20.46 |
| 70               | SAR  | 143.60  | 18.58 | 166.93  | 19.27 | 150.52  | 18.93 | 159.28  | 19.25 | 154.28  | 19.17 |
| 105-106          | Dose | 87.15   | 17.99 | 97.51   | 18.47 | 91.32   | 18.18 | 99.86   | 18.53 | 101.10  | 19.12 |
| 129              | SAR  | 42.03   | 19.87 | 39.31   | 20.28 | 42.43   | 18.17 | 34.56   | 20.05 | 40.71   | 20.43 |
| 130              | SAR  | 138.57  | 17.76 | 147.99  | 18.68 | 142.98  | 19.29 | 150.27  | 18.49 | 148.09  | 20.41 |
| 165-166          | Dose | 89.15   | 18.04 | 86.24   | 18.93 | 92.27   | 19.51 | 88.26   | 18.65 | 96.64   | 18.59 |
| 187              | SAR  | 36.33   | 20.01 | 37.75   | 20.06 | 51.64   | 20.84 | 33.77   | 19.84 | 44.69   | 21.23 |
| 188              | SAR  | 138.95  | 18.61 | 148.22  | 19.36 | 130.27  | 19.40 | 144.81  | 20.04 | 142.98  | 20.05 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 6 | Can 7 | Can 8 | Can 9 | Can 10 |
|-----------|-------|-------|-------|-------|--------|
| 16-20     | 706   | 699   | 689   | 705   | 674    |
| 199-203   | 776   | 767   | 778   | 788   | 761    |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7503

Theoretical Water Content: 1000ppm

Analyst: Brook White

| Actuation Number | Test | Inverted            |                    |         |       |         |       |         |       |         |       |
|------------------|------|---------------------|--------------------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1               |                    | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
| 25               | SAR  | Alb SO4             | IPB                | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 26               | SAR  | 119.53              | 21.41              | 115.26  | 20.43 | 120.54  | 21.28 | 119.53  | 21.11 | 119.22  | 20.97 |
| 33               | SAR  | 117.87              | 20.41              | 118.87  | 21.03 | 120.08  | 20.94 | 127.94  | 22.08 | 122.63  | 21.68 |
| 34               | SAR  | 66.45               | 20.06              | 61.57   | 20.56 | 68.92   | 20.67 | 75.56   | 21.64 | 64.60   | 20.49 |
| 45-46            | Dose | 154.57              | 21.19              | 149.77  | 20.11 | 147.67  | 20.40 | 149.92  | 20.42 | 157.16  | 20.72 |
| 69               | SAR  | 103.75              | 20.11              | 102.91  | 19.73 | 106.51  | 20.53 | 111.15  | 20.83 | 110.27  | 20.56 |
| 70               | SAR  | 60.79               | 19.61              | 55.17   | 19.58 | 60.55   | 19.83 | 62.48   | 20.96 | 67.43   | 19.78 |
| 105-106          | Dose | 152.72              | 20.44              | 144.75  | 19.98 | 156.18  | 20.81 | 157.57  | 20.73 | 144.23  | 20.33 |
| 129              | SAR  | 112.72              | 20.10              | 110.56  | 20.10 | 113.97  | 20.68 | 114.57  | 20.73 | 113.18  | 20.34 |
| 130              | SAR  | 72.77 <sup>1</sup>  | 20.19 <sup>1</sup> | 57.44   | 19.53 | 63.86   | 19.95 | 65.41   | 20.83 | 66.83   | 19.98 |
| 165-166          | Dose | 218.36 <sup>2</sup> | 51.20 <sup>2</sup> | 141.27  | 20.17 | 143.49  | 20.80 | 155.02  | 22.09 | 147.87  | 20.91 |
| 187              | SAR  | 103.11              | 19.98              | 101.52  | 19.97 | 108.29  | 20.80 | 109.34  | 20.98 | 103.87  | 20.32 |
| 188              | SAR  | 89.85               | 19.85              | 63.00   | 19.70 | 67.44   | 20.67 | 78.93   | 22.96 | 70.80   | 20.65 |
|                  |      | 135.42              | 21.69              | 138.18  | 20.71 | 140.29  | 21.73 | 143.79  | 21.50 | 145.14  | 22.12 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 892   | 911   | 919   | 918   | 880   |
| 199-203   | 1019  | 995   | 1018  | 1033  | 983   |

<sup>1</sup>Data was confirmed by reinjection due to out of trend results for the second actuation of this test. Results of 71.96 and 19.96 µg obtained for Albuterol Sulfate and Ipratropium Bromide respectively.

<sup>2</sup>Data was confirmed by reinjection since it was out of trend with the data set. Results of 214.59 and 49.91 µg were obtained for Albuterol Sulfate and Ipratropium Bromide respectively. Analytical error is suspected but could not be confirmed. The sponsor requested to report the data with no further investigation.



# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7503

Theoretical Water Content: 1000ppm

Analyst: Nichole Byrum

| Actuation Number | Test | Upright |       |         |       |         |                   |         |       |         |       |
|------------------|------|---------|-------|---------|-------|---------|-------------------|---------|-------|---------|-------|
|                  |      | Can 6   |       | Can 7   |       | Can 8   |                   | Can 9   |       | Can 10  |       |
| 25               | SAR  | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB               | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 26               | SAR  | 119.34  | 21.03 | 117.10  | 21.04 | 117.19  | 20.91             | 116.11  | 20.80 | 115.67  | 20.69 |
| 33               | SAR  | 123.22  | 21.36 | 116.02  | 20.31 | 117.45  | 20.44             | 117.78  | 20.58 | 118.97  | 20.92 |
| 34               | SAR  | 108.61  | 20.57 | 114.23  | 20.03 | 113.92  | 19.80             | 114.38  | 20.00 | 110.03  | 19.08 |
| 45-46            | Dose | 132.58  | 20.98 | 117.31  | 19.03 | 114.35  | 19.27             | 121.27  | 19.69 | 120.60  | 20.04 |
| 69               | SAR  | 102.75  | 20.12 | 106.17  | 19.35 | 110.84  | 18.77             | 109.84  | 19.40 | 106.21  | 18.16 |
| 70               | SAR  | 39.01   | 19.73 | 42.26   | 19.64 | 43.31   | 19.58             | 49.47   | 20.20 | 39.69   | 19.38 |
| 105-106          | Dose | 168.30  | 21.82 | 160.79  | 19.54 | 44.79   | 5.41 <sup>1</sup> | 152.80  | 19.17 | 148.38  | 18.04 |
| 129              | SAR  | 116.50  | 20.06 | 114.20  | 19.35 | 108.99  | 19.00             | 113.06  | 19.37 | 110.08  | 18.83 |
| 130              | SAR  | 56.87   | 21.02 | 96.09   | 19.72 | 107.59  | 20.20             | 105.46  | 20.43 | 87.33   | 20.02 |
| 165-166          | Dose | 152.62  | 20.61 | 128.39  | 19.42 | 117.70  | 19.45             | 125.52  | 19.90 | 134.88  | 18.48 |
| 187              | SAR  | 108.03  | 20.49 | 117.27  | 19.70 | 105.54  | 19.52             | 111.27  | 19.09 | 111.82  | 18.53 |
| 188              | SAR  | 88.56   | 20.95 | 42.13   | 20.80 | 53.94   | 20.73             | 103.14  | 20.27 | 42.25   | 19.97 |
|                  |      | 128.73  | 20.91 | 147.33  | 19.53 | 146.23  | 19.64             | 117.17  | 20.09 | 147.48  | 20.06 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 6 | Can 7 | Can 8 | Can 9 | Can 10 |
|-----------|-------|-------|-------|-------|--------|
| 16-20     | 936   | 909   | 886   | 915   | 900    |
| 199-203   | 1045  | 1018  | 1013  | 1020  | 998    |

<sup>1</sup> This value was lower than expected. However, there was no follow-up investigation of this result. This result was presented to the sponsor. The sponsor had no additional comments.

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7504

Theoretical Water Content: 1200ppm

Analyst: Dan Milka

| Actuation Number | Test | Inverted |       |         |       |         |       |         |       |         |       |
|------------------|------|----------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1    |       | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
| 25               | SAR  | Alb SO4  | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 26               | SAR  | 116.12   | 20.93 | 116.98  | 20.70 | 117.61  | 20.87 | 116.53  | 20.71 | 115.57  | 20.60 |
| 33               | SAR  | 116.88   | 20.19 | 117.37  | 20.54 | 117.99  | 20.42 | 117.06  | 20.56 | 117.26  | 20.44 |
| 34               | SAR  | 115.47   | 20.64 | 120.17  | 20.28 | 117.21  | 20.08 | 116.66  | 20.77 | 117.19  | 23.65 |
| 45-46            | Dose | 117.83   | 20.07 | 114.54  | 20.87 | 118.94  | 20.22 | 120.28  | 20.56 | 115.64  | 20.19 |
| 69               | SAR  | 111.40   | 19.73 | 112.34  | 26.72 | 109.13  | 19.19 | 113.77  | 20.06 | 114.25  | 19.97 |
| 70               | SAR  | 85.73    | 20.12 | 100.71  | 19.82 | 86.07   | 20.13 | 88.78   | 21.01 | 94.15   | 19.54 |
| 105-106          | Dose | 138.10   | 20.16 | 123.70  | 20.33 | 139.96  | 20.57 | 131.06  | 19.38 | 132.44  | 20.56 |
| 129              | SAR  | 110.05   | 19.45 | 110.63  | 19.63 | 112.42  | 19.86 | 110.90  | 19.74 | 110.94  | 19.75 |
| 130              | SAR  | 81.29    | 19.86 | 80.29   | 19.96 | 88.14   | 20.50 | 82.38   | 19.87 | 88.28   | 20.43 |
| 165-166          | Dose | 136.33   | 19.87 | 138.94  | 20.34 | 140.05  | 20.22 | 139.48  | 20.73 | 145.92  | 21.22 |
| 187              | SAR  | 107.59   | 19.89 | 107.69  | 19.87 | 108.52  | 19.97 | 108.39  | 20.14 | 110.74  | 20.19 |
| 188              | SAR  | 113.81   | 20.29 | 118.65  | 20.79 | 102.22  | 21.08 | 97.04   | 20.54 | 101.05  | 20.69 |
|                  |      | 117.03   | 20.89 | 111.13  | 20.59 | 120.00  | 20.58 | 127.32  | 21.03 | 117.76  | 20.01 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 1081  | 1102  | 1066  | 1064  | 1000  |
| 199-203   | 1168  | 1169  | 1155  | 1160  | 1099  |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7504

Theoretical Water Content: 1200ppm

Analyst: Kaitryn Johnson

| Activation Number | Test | Upright |       |         |       |         |       |         |       |         |        |
|-------------------|------|---------|-------|---------|-------|---------|-------|---------|-------|---------|--------|
|                   |      | Can 6   | Can 6 | Can 7   | Can 7 | Can 8   | Can 8 | Can 9   | Can 9 | Can 10  | Can 10 |
| 25                | SAR  | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB    |
| 26                | SAR  | 116.10  | 21.10 | 114.03  | 20.41 | 112.53  | 20.47 | 114.66  | 20.89 | 116.74  | 20.97  |
| 33                | SAR  | 113.22  | 19.97 | 116.02  | 20.73 | 115.40  | 20.57 | 118.59  | 21.44 | 113.78  | 20.05  |
| 34                | SAR  | 48.49   | 20.05 | 67.76   | 19.86 | 43.53   | 18.90 | 55.28   | 19.82 | 58.49   | 20.20  |
| 45-46             | Dose | 158.76  | 18.40 | 155.17  | 18.98 | 170.51  | 19.45 | 162.81  | 18.78 | 157.29  | 18.99  |
| 69                | SAR  | 100.66  | 17.79 | 104.98  | 18.43 | 98.69   | 17.57 | 107.29  | 19.38 | 103.50  | 19.15  |
| 70                | SAR  | 40.69   | 19.14 | 40.95   | 19.10 | 50.42   | 18.85 | 49.55   | 19.44 | 41.17   | 18.60  |
| 105-106           | Dose | 152.82  | 17.78 | 155.40  | 17.44 | 147.61  | 17.45 | 155.81  | 19.48 | 167.65  | 18.68  |
| 129               | SAR  | 103.29  | 17.92 | 103.79  | 18.25 | 102.22  | 17.79 | 105.39  | 18.17 | 102.98  | 18.32  |
| 130               | SAR  | 52.51   | 19.10 | 59.34   | 19.88 | 57.51   | 19.25 | 75.05   | 19.87 | 51.54   | 19.51  |
| 165-166           | Dose | 138.19  | 17.26 | 144.92  | 19.03 | 138.06  | 17.93 | 133.71  | 19.07 | 145.81  | 18.57  |
| 187               | SAR  | 93.54   | 16.06 | 96.28   | 16.46 | 93.60   | 17.05 | 99.05   | 17.39 | 99.44   | 17.26  |
| 188               | SAR  | 31.05   | 19.40 | 34.49   | 20.21 | 37.58   | 18.60 | 43.68   | 20.52 | 31.99   | 20.24  |
|                   |      | 147.05  | 18.96 | 151.85  | 19.20 | 144.52  | 18.83 | 146.22  | 19.65 | 148.92  | 19.78  |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Activation Reproducibility

## Water Content (ppm)

| Activation | Can 6 | Can 7 | Can 8 | Can 9 | Can 10 |
|------------|-------|-------|-------|-------|--------|
| 16-20      | 1063  | 1066  | 1060  | 1071  | 1082   |
| 199-203    | 1159  | 1173  | 1154  | 1175  | 1184   |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7505

Theoretical Water Content: 1500ppm

Analyst: Brook White

| Actuation Number | Test | Inverted |       |         |       |         |       |         |       |         |       |
|------------------|------|----------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1    |       | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
|                  |      | Alb SO4  | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 116.21   | 20.67 | 115.98  | 20.74 | 119.88  | 21.38 | 126.79  | 22.26 | 118.66  | 21.10 |
| 26               | SAR  | 118.46   | 20.67 | 116.09  | 20.10 | 121.20  | 21.20 | 122.00  | 21.34 | 116.72  | 20.73 |
| 33               | SAR  | 112.55   | 19.84 | 116.28  | 20.82 | 114.67  | 20.34 | 113.06  | 20.31 | 114.61  | 20.86 |
| 34               | SAR  | 117.06   | 20.50 | 111.72  | 19.32 | 118.41  | 20.54 | 110.90  | 19.19 | 115.39  | 20.39 |
| 45-46            | Dose | 110.45   | 19.80 | 111.29  | 19.78 | 114.27  | 20.55 | 112.32  | 19.91 | 113.49  | 20.17 |
| 69               | SAR  | 99.24    | 19.26 | 108.44  | 20.24 | 99.74   | 20.02 | 105.97  | 20.25 | 101.45  | 20.17 |
| 70               | SAR  | 122.24   | 20.06 | 119.12  | 20.08 | 124.01  | 20.35 | 118.92  | 19.94 | 124.88  | 20.31 |
| 105-106          | Dose | 111.73   | 20.21 | 113.12  | 20.06 | 113.20  | 20.48 | 112.03  | 19.96 | 114.97  | 20.84 |
| 129              | SAR  | 113.98   | 20.14 | 113.93  | 19.94 | 112.79  | 20.78 | 110.89  | 20.07 | 108.77  | 19.69 |
| 130              | SAR  | 118.36   | 21.38 | 115.40  | 20.20 | 116.83  | 20.81 | 115.00  | 20.45 | 125.57  | 22.64 |
| 165-166          | Dose | 107.47   | 19.95 | 112.14  | 20.24 | 110.44  | 20.25 | 109.26  | 20.24 | 109.89  | 20.29 |
| 187              | SAR  | 117.19   | 20.72 | 119.51  | 20.61 | 110.54  | 20.51 | 107.43  | 20.37 | 119.79  | 20.74 |
| 188              | SAR  | 107.92   | 20.35 | 114.33  | 20.87 | 122.09  | 22.64 | 112.39  | 21.22 | 117.01  | 22.25 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 1381  | 1430  | 1361  | 1413  | 1380  |
| 199-203   | 1488  | 1488  | 1438  | 1503  | 1486  |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7505

Theoretical Water Content: 1500ppm

Analyst: Nichole Bynum

| Actuation Number | Test | Upright |       |         |       |         |       |         |       |         |       |
|------------------|------|---------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 6   |       | Can 7   |       | Can 8   |       | Can 9   |       | Can 10  |       |
|                  |      | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 116.73  | 20.84 | 117.75  | 20.66 | 119.88  | 20.65 | 117.66  | 20.83 | 125.01  | 21.96 |
| 26               | SAR  | 113.69  | 20.06 | 115.37  | 20.18 | 118.13  | 20.57 | 125.89  | 21.77 | 124.74  | 21.81 |
| 33               | SAR  | 119.64  | 22.47 | 116.18  | 19.72 | 118.66  | 20.09 | 119.61  | 20.44 | 116.61  | 20.40 |
| 34               | SAR  | 119.72  | 19.94 | 117.46  | 19.94 | 118.13  | 19.56 | 119.28  | 19.61 | 121.44  | 20.13 |
| 45-46            | Dose | 114.78  | 19.30 | 114.70  | 19.16 | 108.98  | 18.65 | 114.77  | 19.74 | 112.26  | 18.47 |
| 69               | SAR  | 112.18  | 20.14 | 116.05  | 19.46 | 111.63  | 19.43 | 114.12  | 20.20 | 114.19  | 19.60 |
| 70               | SAR  | 122.24  | 20.85 | 117.55  | 19.33 | 110.99  | 18.25 | 118.92  | 20.47 | 118.45  | 18.81 |
| 105-106          | Dose | 116.16  | 19.95 | 114.83  | 19.56 | 112.33  | 19.05 | 112.30  | 20.42 | 114.25  | 19.85 |
| 129              | SAR  | 113.66  | 20.70 | 109.67  | 19.65 | 122.59  | 19.86 | 116.34  | 20.45 | 112.39  | 19.84 |
| 130              | SAR  | 113.45  | 20.18 | 121.31  | 19.92 | 85.89   | 14.92 | 112.48  | 19.30 | 121.42  | 20.02 |
| 165-166          | Dose | 109.10  | 20.13 | 108.94  | 19.64 | 114.02  | 18.91 | 110.36  | 20.11 | 115.55  | 17.60 |
| 187              | SAR  | 97.28   | 20.74 | 99.90   | 19.93 | 105.47  | 20.31 | 107.48  | 20.96 | 104.49  | 20.53 |
| 188              | SAR  | 126.94  | 21.01 | 116.29  | 19.90 | 111.20  | 19.25 | 111.98  | 20.10 | 110.54  | 20.38 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 6 | Can 7 | Can 8 | Can 9 | Can 10 |
|-----------|-------|-------|-------|-------|--------|
| 16-20     | 1367  | 1393  | 1369  | 1428  | 1381   |
| 199-203   | 1491  | 1492  | 1471  | 1562  | 1498   |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7506

Theoretical Water Content: 2500ppm

Analyst: Dan Milka

| Actuation Number | Test | Inverted |       |         |       |         |       |         |       |         |       |
|------------------|------|----------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1    |       | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
| 25               | SAR  | Alb SO4  | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 26               | SAR  | 122.31   | 20.76 | 124.15  | 21.38 | 119.40  | 20.64 | 119.42  | 20.36 | 119.08  | 20.20 |
| 33               | SAR  | 113.19   | 20.34 | 116.23  | 20.70 | 116.41  | 20.91 | 112.91  | 20.41 | 114.47  | 20.75 |
| 34               | SAR  | 119.83   | 20.11 | 122.32  | 21.04 | 118.66  | 20.20 | 117.66  | 20.03 | 116.36  | 20.43 |
| 45-46            | Dose | 111.77   | 20.44 | 115.00  | 20.53 | 116.50  | 21.09 | 111.81  | 20.13 | 111.95  | 20.65 |
| 69               | SAR  | 113.79   | 19.84 | 115.73  | 20.62 | 114.49  | 20.17 | 110.62  | 19.64 | 111.29  | 20.07 |
| 70               | SAR  | 116.56   | 19.96 | 122.37  | 20.74 | 120.39  | 20.43 | 117.88  | 20.32 | 120.73  | 20.49 |
| 105-106          | Dose | 114.39   | 20.46 | 112.99  | 20.27 | 112.45  | 20.60 | 111.14  | 19.64 | 108.17  | 20.23 |
| 129              | SAR  | 115.43   | 20.35 | 117.74  | 20.53 | 115.95  | 20.45 | 114.04  | 19.91 | 112.84  | 19.81 |
| 130              | SAR  | 118.77   | 20.26 | 119.46  | 20.71 | 118.05  | 20.49 | 116.71  | 19.78 | 118.50  | 20.18 |
| 165-166          | Dose | 117.04   | 21.16 | 113.68  | 21.26 | 115.08  | 21.15 | 108.10  | 20.10 | 108.95  | 20.17 |
| 187              | SAR  | 109.04   | 20.06 | 111.09  | 20.55 | 108.86  | 19.84 | 107.36  | 19.43 | 106.84  | 19.87 |
| 188              | SAR  | 112.35   | 20.52 | 114.37  | 20.90 | 113.56  | 20.83 | 112.29  | 20.10 | 113.38  | 20.71 |
|                  |      | 109.54   | 20.86 | 112.00  | 21.27 | 110.88  | 21.09 | 106.57  | 20.69 | 106.19  | 21.48 |

Alb SO4 = Alburexol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 2368  | 2372  | 2336  | 2313  | 2370  |
| 199-203   | 2502  | 2483  | 2466  | 2479  | 2465  |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7506

Theoretical Water Content: 2500ppm

Analyst: Kathryn Johnson

| Actuation Number | Test | Upright |       |         |       |         |       |         |       |         |       |
|------------------|------|---------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 6   |       | Can 7   |       | Can 8   |       | Can 9   |       | Can 10  |       |
|                  |      | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 118.31  | 20.22 | 121.52  | 20.71 | 122.54  | 21.01 | 120.44  | 21.14 | 119.51  | 20.58 |
| 26               | SAR  | 116.30  | 21.63 | 115.59  | 21.24 | 115.78  | 20.96 | 115.26  | 21.33 | 118.52  | 21.78 |
| 33               | SAR  | 118.15  | 19.67 | 123.20  | 20.24 | 108.74  | 18.68 | 112.74  | 18.66 | 116.39  | 19.08 |
| 34               | SAR  | 110.16  | 19.68 | 111.55  | 19.23 | 123.30  | 19.05 | 118.92  | 19.67 | 112.31  | 18.76 |
| 45-46            | Dose | 108.85  | 19.69 | 115.19  | 20.10 | 109.42  | 18.35 | 110.01  | 18.79 | 112.51  | 19.01 |
| 69               | SAR  | 114.10  | 19.92 | 116.72  | 20.56 | 109.25  | 19.37 | 110.22  | 18.85 | 112.76  | 20.10 |
| 70               | SAR  | 108.20  | 18.66 | 109.74  | 18.78 | 108.62  | 17.54 | 110.10  | 18.58 | 106.14  | 19.68 |
| 105-106          | Dose | 103.46  | 18.53 | 106.70  | 19.25 | 103.64  | 18.18 | 104.61  | 18.86 | 102.27  | 18.19 |
| 129              | SAR  | 105.51  | 19.90 | 106.09  | 19.84 | 104.55  | 19.93 | 105.78  | 19.51 | 107.21  | 19.65 |
| 130              | SAR  | 104.97  | 19.09 | 106.68  | 19.73 | 111.66  | 18.76 | 106.94  | 19.75 | 108.97  | 19.96 |
| 165-166          | Dose | 95.47   | 18.82 | 99.41   | 19.17 | 96.08   | 17.84 | 98.80   | 19.22 | 94.20   | 18.74 |
| 187              | SAR  | 100.26  | 21.00 | 101.76  | 20.59 | 97.34   | 20.18 | 98.89   | 20.03 | 93.45   | 20.01 |
| 188              | SAR  | 95.69   | 20.01 | 95.40   | 20.51 | 102.57  | 18.92 | 97.96   | 20.08 | 92.92   | 19.25 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 6 | Can 7 | Can 8 | Can 9 | Can 10 |
|-----------|-------|-------|-------|-------|--------|
| 16-20     | 2362  | 2356  | 2327  | 2306  | 2346   |
| 199-203   | 2473  | 2499  | 2436  | 2460  | 2488   |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7507

Theoretical Water Content: 3000ppm

Analyst: Brook White

| Actuation Number | Test | Inverted |       |         |       |         |       |         |       |         |       |
|------------------|------|----------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1    |       | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
|                  |      | Alb SO4  | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 116.79   | 20.63 | 118.78  | 20.84 | 114.65  | 19.88 | 115.17  | 20.26 | 114.81  | 20.36 |
| 26               | SAR  | 116.64   | 20.34 | 114.57  | 20.24 | 113.22  | 20.02 | 117.60  | 21.12 | 114.75  | 20.95 |
| 33               | SAR  | 116.12   | 20.18 | 114.51  | 19.30 | 111.95  | 19.23 | 116.47  | 20.44 | 114.74  | 20.25 |
| 34               | SAR  | 115.87   | 20.62 | 113.18  | 19.59 | 109.45  | 19.21 | 113.72  | 20.14 | 112.13  | 19.84 |
| 45-46            | Dose | 113.55   | 20.27 | 111.76  | 19.86 | 110.04  | 19.66 | 112.80  | 20.11 | 110.32  | 19.90 |
| 69               | SAR  | 109.91   | 19.33 | 114.65  | 20.07 | 109.43  | 18.77 | 114.66  | 19.86 | 113.59  | 20.15 |
| 70               | SAR  | 110.79   | 20.15 | 111.36  | 19.66 | 108.25  | 19.50 | 108.31  | 19.41 | 112.30  | 19.95 |
| 105-106          | Dose | 106.24   | 19.04 | 105.50  | 18.97 | 104.43  | 18.94 | 108.74  | 19.63 | 105.61  | 19.64 |
| 129              | SAR  | 111.53   | 19.66 | 112.73  | 20.08 | 108.28  | 19.39 | 112.30  | 19.79 | 111.44  | 19.67 |
| 130              | SAR  | 111.51   | 20.43 | 112.78  | 20.33 | 107.66  | 19.59 | 110.29  | 20.15 | 111.13  | 20.54 |
| 165-166          | Dose | 104.08   | 19.03 | 105.03  | 19.81 | 102.79  | 19.27 | 106.09  | 20.13 | 101.12  | 19.55 |
| 187              | SAR  | 105.30   | 19.64 | 106.67  | 20.02 | 105.85  | 19.30 | 108.70  | 20.34 | 104.53  | 20.37 |
| 188              | SAR  | 108.80   | 21.15 | 105.98  | 20.34 | 102.38  | 20.16 | 109.39  | 21.92 | 106.77  | 21.39 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 2880  | 2897  | 2847  | 2812  | 3071  |
| 199-203   | 3072  | 3059  | 3016  | 3040  | 3304  |



# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7507

Theoretical Water Content: 3000ppm

Analyst: Kathryn Johnson

| Actuation Number | Test | Upright |       |         |       |                     |                    |         |       |         |       |
|------------------|------|---------|-------|---------|-------|---------------------|--------------------|---------|-------|---------|-------|
|                  |      | Can 6   |       | Can 7   |       | Can 8               |                    | Can 9   |       | Can 10  |       |
| 25               | SAR  | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4             | IPB                | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 26               | SAR  | 115.63  | 20.16 | 113.86  | 19.80 | 116.54              | 20.24              | 112.98  | 19.99 | 115.98  | 20.94 |
| 33               | SAR  | 107.23  | 19.69 | 107.73  | 19.16 | 115.93              | 20.91              | 111.69  | 20.12 | 111.83  | 20.05 |
| 34               | SAR  | 115.88  | 19.98 | 114.25  | 19.05 | 114.43              | 19.24              | 113.75  | 18.92 | 114.69  | 19.39 |
| 45-46            | Dose | 105.62  | 17.80 | 108.54  | 17.95 | 114.01              | 18.69              | 112.02  | 18.30 | 111.47  | 18.97 |
| 69               | SAR  | 103.31  | 17.63 | 104.60  | 17.76 | 109.46              | 18.27              | 105.77  | 17.84 | 104.59  | 17.82 |
| 70               | SAR  | 101.08  | 19.01 | 106.33  | 18.83 | 107.76              | 19.42              | 108.06  | 19.17 | 107.30  | 19.12 |
| 105-106          | Dose | 108.12  | 17.37 | 106.76  | 17.85 | 113.24 <sup>1</sup> | 18.66 <sup>1</sup> | 108.54  | 19.06 | 104.75  | 18.35 |
| 129              | SAR  | 101.02  | 17.35 | 102.02  | 18.07 | 103.96              | 17.68              | 102.51  | 18.06 | 101.53  | 18.08 |
| 130              | SAR  | 103.08  | 19.77 | 103.14  | 18.08 | 102.88              | 19.45              | 108.60  | 19.57 | 104.96  | 19.51 |
| 165-166          | Dose | 106.76  | 17.27 | 108.04  | 18.54 | 115.61              | 19.45              | 102.71  | 18.65 | 103.97  | 18.55 |
| 187              | SAR  | 90.88   | 12.85 | 70.26   | 8.52  | 87.87               | 10.13              | 89.64   | 12.38 | 89.04   | 13.95 |
| 188              | SAR  | 96.45   | 19.42 | 90.68   | 19.05 | 87.57               | 19.62              | 87.02   | 19.59 | 86.35   | 19.42 |
|                  |      | 88.82   | 18.99 | 90.26   | 18.84 | 100.53              | 20.92              | 91.29   | 19.42 | 80.36   | 19.44 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 6 | Can 7 | Can 8 | Can 9 | Can 10 |
|-----------|-------|-------|-------|-------|--------|
| 16-20     | 2795  | 2855  | 2831  | 2805  | 2841   |
| 199-203   | 2956  | 3049  | 3004  | 2976  | 3008   |

<sup>1</sup> Analyst observed that this actuation sounded weaker than normal.

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7508

Theoretical Water Content: 3500ppm

Analyst: Brook White

| Actuation Number | Test | Inverted |       |         |       |         |       |         |       |         |       |
|------------------|------|----------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1    |       | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
|                  |      | Alb SO4  | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 119.26   | 21.31 | 114.26  | 20.00 | 113.72  | 20.20 | 109.36  | 19.76 | 111.63  | 19.84 |
| 26               | SAR  | 116.90   | 20.55 | 104.49  | 19.16 | 108.69  | 19.55 | 115.32  | 20.69 | 112.33  | 19.77 |
| 33               | SAR  | 109.74   | 19.54 | 115.38  | 20.34 | 110.26  | 19.76 | 111.04  | 19.22 | 114.81  | 20.06 |
| 34               | SAR  | 114.78   | 20.54 | 109.11  | 20.00 | 110.68  | 19.26 | 113.07  | 20.00 | 114.46  | 20.74 |
| 45-46            | Dose | 108.77   | 19.61 | 106.58  | 19.62 | 106.91  | 19.52 | 108.66  | 20.02 | 108.71  | 19.72 |
| 69               | SAR  | 110.32   | 19.31 | 111.13  | 20.05 | 103.45  | 19.29 | 108.26  | 19.49 | 110.96  | 19.65 |
| 70               | SAR  | 110.29   | 20.28 | 109.88  | 19.94 | 108.21  | 20.05 | 108.26  | 20.34 | 113.82  | 20.50 |
| 105-106          | Dose | 102.54   | 18.78 | 105.92  | 19.36 | 103.81  | 19.46 | 104.38  | 19.32 | 105.34  | 19.46 |
| 129              | SAR  | 112.03   | 20.40 | 108.01  | 19.88 | 107.55  | 20.08 | 105.69  | 19.45 | 109.58  | 20.01 |
| 130              | SAR  | 110.87   | 20.93 | 109.67  | 20.58 | 108.77  | 20.84 | 108.23  | 19.97 | 106.11  | 19.70 |
| 165-166          | Dose | 101.92   | 19.77 | 98.75   | 19.31 | 97.24   | 18.46 | 100.21  | 19.81 | 100.89  | 19.01 |
| 187              | SAR  | 102.85   | 20.05 | 104.49  | 20.09 | 106.31  | 20.19 | 99.75   | 19.07 | 105.50  | 19.50 |
| 188              | SAR  | 104.69   | 21.04 | 99.83   | 20.07 | 103.75  | 21.25 | 101.51  | 21.24 | 99.90   | 20.54 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuations | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|------------|-------|-------|-------|-------|-------|
| 16-20      | 3416  | 3352  | 3367  | 3320  | 3302  |
| 199-203    | 3533  | 3556  | 3558  | 3540  | 3551  |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7508

Theoretical Water Content: 3500ppm

Analyst: Kathryn Johnson

| Actuation Number | Test | Upright |       |         |       |         |       |                    |       |         |       |
|------------------|------|---------|-------|---------|-------|---------|-------|--------------------|-------|---------|-------|
|                  |      | Can 6   |       | Can 7   |       | Can 8   |       | Can 9 <sup>1</sup> |       | Can 10  |       |
| 25               | SAR  | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4            | IPB   | Alb SO4 | IPB   |
| 26               | SAR  | 112.02  | 20.13 | 114.41  | 19.95 | 115.22  | 19.99 | 117.50             | 20.48 | 116.70  | 20.62 |
| 33               | SAR  | 108.44  | 19.80 | 114.58  | 21.22 | 109.67  | 20.27 | 109.60             | 20.23 | 107.58  | 19.57 |
| 34               | SAR  | 109.10  | 18.51 | 109.33  | 18.95 | 117.84  | 20.51 | 112.49             | 18.89 | 110.15  | 19.19 |
| 45-46            | Dose | 110.44  | 18.09 | 100.80  | 19.74 | 106.24  | 19.27 | 111.46             | 18.80 | 116.55  | 18.97 |
| 69               | SAR  | 105.40  | 18.21 | 99.55   | 17.22 | 109.68  | 19.51 | 106.55             | 17.97 | 108.31  | 18.30 |
| 70               | SAR  | 101.42  | 18.98 | 102.63  | 18.04 | 104.06  | 18.76 | 104.88             | 17.15 | 103.54  | 19.23 |
| 105-106          | Dose | 105.25  | 17.95 | 103.34  | 17.65 | 114.88  | 18.83 | 105.91             | 17.79 | 117.73  | 19.08 |
| 129              | SAR  | 101.56  | 18.03 | 101.28  | 17.68 | 105.87  | 18.39 | 103.37             | 17.85 | 104.14  | 17.96 |
| 130              | SAR  | 96.51   | 19.09 | 98.10   | 17.55 | 101.57  | 18.96 | 99.14              | 18.12 | 103.08  | 19.85 |
| 165-166          | Dose | 106.27  | 18.53 | 98.85   | 17.92 | 112.53  | 19.76 | 101.08             | 19.14 | 107.76  | 18.90 |
| 187              | SAR  | 88.29   | 12.23 | 85.45   | 11.54 | 92.06   | 13.77 | 91.10              | 13.67 | 92.48   | 14.08 |
| 188              | SAR  | 90.87   | 19.18 | 94.74   | 19.23 | 98.05   | 19.83 | 87.74              | 19.18 | 91.15   | 20.00 |
|                  |      | 84.44   | 19.42 | 85.29   | 19.97 | 86.82   | 19.65 | 81.18              | 19.60 | 83.01   | 19.58 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuations | Can 6 | Can 7 | Can 8 | Can 9 <sup>1</sup> | Can 10 |
|------------|-------|-------|-------|--------------------|--------|
| 16-20      | 3360  | 3381  | 3404  | 3396               | 3365   |
| 199-203    | 3528  | 3360  | 3572  | 3567               | 3551   |

<sup>1</sup>During initial water content testing, 6 actuations may have been made (16-21 actuations), Sponsor was notified.

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7509

Theoretical Water Content: 4000ppm

Analyst: Brook White

| Actuation Number | Test | Inverted |       |         |       |         |       |         |       |         |       |
|------------------|------|----------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1    |       | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
| 25               | SAR  | Alb SO4  | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 26               | SAR  | 111.00   | 20.05 | 111.12  | 20.08 | 107.69  | 19.85 | 109.37  | 19.78 | 109.83  | 20.26 |
| 33               | SAR  | 114.93   | 19.82 | 108.69  | 19.84 | 105.13  | 19.80 | 111.82  | 20.29 | 109.57  | 20.61 |
| 34               | SAR  | 111.64   | 19.64 | 108.50  | 19.29 | 107.54  | 19.66 | 109.40  | 19.63 | 112.92  | 20.30 |
| 45-46            | Dose | 108.96   | 19.71 | 103.73  | 19.06 | 103.49  | 19.75 | 110.56  | 19.85 | 109.26  | 20.35 |
| 69               | SAR  | 105.70   | 19.34 | 107.09  | 19.33 | 103.04  | 18.84 | 107.71  | 19.52 | 108.75  | 20.06 |
| 70               | SAR  | 107.67   | 19.19 | 108.24  | 19.42 | 103.84  | 19.33 | 103.61  | 19.28 | 109.93  | 19.81 |
| 105-106          | Dose | 110.99   | 19.80 | 106.52  | 20.27 | 100.93  | 18.86 | 109.04  | 20.61 | 107.58  | 20.24 |
| 129              | SAR  | 103.56   | 18.83 | 102.73  | 19.31 | 101.74  | 19.22 | 104.33  | 19.37 | 107.92  | 19.91 |
| 130              | SAR  | 108.69   | 19.43 | 106.55  | 19.37 | 103.40  | 19.04 | 106.54  | 20.05 | 110.48  | 20.40 |
| 165-166          | Dose | 105.58   | 20.83 | 105.30  | 20.37 | 97.61   | 18.90 | 102.24  | 19.86 | 107.98  | 20.88 |
| 187              | SAR  | 93.00    | 18.04 | 95.39   | 19.28 | 91.81   | 18.85 | 95.40   | 19.48 | 99.05   | 19.48 |
| 188              | SAR  | 99.14    | 20.18 | 96.00   | 19.67 | 97.45   | 19.73 | 98.47   | 19.44 | 97.47   | 19.54 |
|                  |      | 94.35    | 20.32 | 93.34   | 20.10 | 95.22   | 20.11 | 96.11   | 21.03 | 99.88   | 21.63 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 3805  | 3822  | 3843  | 3754  | 3811  |
| 199-203   | 4045  | 4003  | 4033  | 3938  | 4014  |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7509

Theoretical Water Content: 4000ppm

Analyst: Kathryn Johnson

| Actuation Number | Test | Upright |       |         |       |         |       |         |       |         |       |
|------------------|------|---------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 6   |       | Can 7   |       | Can 8   |       | Can 9   |       | Can 10  |       |
|                  |      | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 116.70  | 20.77 | 110.63  | 19.79 | 109.67  | 19.99 | 108.48  | 19.93 | 113.94  | 20.11 |
| 26               | SAR  | 110.22  | 20.49 | 110.25  | 20.45 | 108.48  | 20.01 | 110.55  | 20.52 | 107.93  | 20.42 |
| 33               | SAR  | 104.40  | 18.46 | 107.39  | 18.49 | 101.89  | 18.89 | 99.49   | 18.78 | 101.71  | 18.83 |
| 34               | SAR  | 112.67  | 17.93 | 112.71  | 17.96 | 113.58  | 17.39 | 110.43  | 17.55 | 109.57  | 17.68 |
| 45-46            | Dose | 107.88  | 18.44 | 104.94  | 17.04 | 103.43  | 17.40 | 102.29  | 17.56 | 101.34  | 17.57 |
| 69               | SAR  | 98.58   | 19.10 | 91.55   | 15.39 | 95.34   | 17.19 | 92.18   | 16.41 | 95.75   | 16.97 |
| 70               | SAR  | 118.85  | 18.97 | 117.62  | 19.28 | 116.26  | 18.73 | 110.24  | 18.37 | 113.12  | 18.90 |
| 105-106          | Dose | 103.70  | 18.06 | 96.35   | 15.50 | 97.83   | 15.05 | 95.71   | 15.62 | 97.11   | 16.02 |
| 129              | SAR  | 99.24   | 19.76 | 93.70   | 18.79 | 94.21   | 18.95 | 88.04   | 18.29 | 94.86   | 19.32 |
| 130              | SAR  | 106.17  | 18.91 | 109.41  | 18.80 | 107.94  | 19.14 | 108.32  | 19.36 | 107.04  | 19.00 |
| 165-166          | Dose | 80.34   | 10.28 | 87.96   | 12.20 | 91.87   | 14.17 | 85.50   | 13.57 | 87.78   | 12.70 |
| 187              | SAR  | 81.95   | 19.26 | 73.03   | 17.19 | 76.95   | 17.66 | 75.68   | 19.52 | 80.47   | 18.37 |
| 188              | SAR  | 74.83   | 19.48 | 85.99   | 18.39 | 87.35   | 19.29 | 66.70   | 18.68 | 89.32   | 20.28 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 6 | Can 7 | Can 8             | Can 9 | Can 10 |
|-----------|-------|-------|-------------------|-------|--------|
| 16-30     | 3815  | 3741  | 3760              | 3829  | 3708   |
| 199-203   | 3995  | 3963  | 3958 <sup>1</sup> | 4030  | 3925   |

<sup>1</sup>This canister was inadvertently tested at actuations 200-204.

## I. Purpose

The purpose of this study is to examine the effects of water content in Ipratropium Bromide/Albuterol Sulfate 0.021/0.120 mg TTV, 10 mL Inhalation Aerosol on testing for Single Actuation and Dose (2 actuations) through life (ca. 200 actuations).

## II. Applicable Regulatory Requirements

This study will be conducted as per GMP (21 CFR 210 and 211) regulations.

## III. Materials

BIPI will manufacture canisters using the same formulation but different water contents. The samples will be grouped as follows.

| AS Testing Number | Lot Number | Target ppm Water |
|-------------------|------------|------------------|
| 02-10-7501        | A1         | Inherent*        |
| 02-10-7502        | A2         | 800*             |
| 02-10-7503        | A3         | 1000*            |
| 02-10-7504        | A4         | 1200*            |
| 02-10-7505        | A5         | 1500*            |
| 02-10-7506        | A6         | 2500*            |
| 02-10-7507        | A7         | 3000             |
| 02-10-7508        | A8         | 3500             |
| 02-10-7509        | B8         | 4000             |

## IV. Study Design

A total of 90 canisters will be tested utilizing the following table. Ten (10) canisters will be selected from each lot, Five (5) to be stored in an upright orientation and Five (5) in an inverted orientation. Canisters will be placed in orientation immediately following initial Water Content testing and will be stored and weighed in their designated orientation throughout the study life.

Water Content testing will be performed at the beginning and end of the study per the following table. Initial Water Content testing will occur no earlier than the Friday prior initiation of the testing schedule. If the initial Water Content is performed on the preceding Friday, then conduct testing as indicated for Beginning Testing Option 1. Otherwise, proceed as indicated in Beginning Testing Option 2.

---

\* Priority sample, test before other samples.

Final Water Content testing may occur as late as the Monday following the end of the study, if necessary. No change to the testing scheme below will be needed if testing is delayed for final Water Content.

For all methods, calculate sample weights as follows:  $W1 - W2 = \text{Sample Weight}$   
W2 is the canister weight recorded immediately following each test point and W1 is the canister weight prior to testing (previous timepoint's W2.)

Wasting actuations will be performed using a jar and baseplate. A delay of no less than 2 hours must occur between completion of previous wasting timepoint and initiation of next wasting timepoint when multiple wasting sessions are required on the same day. The firing and shaking of canisters is to be performed as outlined in the Uniformity of Dose procedure for wasting canisters. A swab will be used to thoroughly clean the interior of the valve stem prior to collecting the final weight for both the Water Content and the Wasting steps. It is critical to clean the valve stem while the interior is still wet from the testing step. Wasting will be performed per the following table.

Samples must equilibrate for 24 hours before any Uniformity of Dose sampling can be performed. The single actuation reproducibility protocol (SAR) requires that two single actuations be collected in separate containers where indicated. Do not shake canisters between actuations for the SAR protocol. Uniformity of Dose samples will be collected per the following table.

Record testing times as follows:

**Water Content testing:** It is not necessary to record testing times.

**Uniformity of Dose/SAR testing:** Record time testing was initiated.

**Wasting:** Record the time testing was initiated and completed.

Beginning Testing Option 1: Water Content testing performed on preceding Friday

| Week 1           |                  |   |
|------------------|------------------|---|
| Day              | Actuation Number | Action  |
| Preceding Friday | 1-5              | Priming   |
|                  | 6-15             | Clearing the Inlet Tube (Water Content)                   |
|                  | 16-20            | Water Content Sample                                      |
| Monday           | 21-24            | Priming   |
|                  | 25 & 26          | 0 hr SAR<br>Shot weight collection                        |
|                  | 27-32            | Actuate to waste after 0 hr SAR<br>Shot weight collection |

Then continue with testing as indicated below for Week 1 Tuesday.

Beginning Testing Option 2: Water Content testing performed on preceding Friday

| Week 1 |                    |   |
|--------|--------------------|---|
| Day    | Actuation Number   | Action  |
| Monday | 1-5 <sup>*</sup>   | Priming   |
|        | 6-15 <sup>*</sup>  | Clearing the Inlet Tube<br>(Water Content)                |
|        | 16-20 <sup>*</sup> | Water Content Sample                                      |
|        | 21 & 22            | 0 hr SAR<br>Shot weight collection                        |
|        | 23-32              | Actuate to waste after 0 hr SAR<br>Shot weight collection |

Then continue with testing as indicated below for Week 1 Tuesday.

| Week 1    |                       |   |
|-----------|-----------------------|---|
| Day       | Actuation Number      | Action  |
| Tuesday   | 33 & 34 <sup>**</sup> | 24 hr SAR<br>Shot weight collection   |
|           | 35-44                 | Actuate to waste after 24 hr SAR<br>Shot weight collection                        |
| Wednesday | 45-46 <sup>**</sup>   | 24 hr Dose<br>Shot weight collection  |
|           | 47-48                 | Actuate to waste after Dose<br>Shot weight collection                             |
|           | 49-52                 | Actuate to waste Mid Day<br>Shot weight collection                                |
|           | 53-56                 | Actuate to waste End Day<br>Shot weight collection                                |
| Thursday  | 57-68                 | Actuate to waste<br>Shot weight collection  |
| Friday    | 69 & 70 <sup>**</sup> | 24 hr SAR<br>Shot weight collection   |
|           | 71-72                 | Actuate to waste after 24 hr SAR<br>(actuators 69 & 70)<br>Shot weight collection |
|           | 73-76                 | Actuate to waste Mid Day<br>Shot weight collection                                |
|           | 77-80                 | Actuate to waste End Day<br>Shot weight collection                                |

<sup>\*</sup> May be performed as early as the Friday prior.

<sup>\*\*</sup> Must occur at least 24 hours from completion of prior wasting timepoint.



| Week 2    |                  |   |
|-----------|------------------|---|
| Day       | Actuation Number | Action  |
| Monday    | 81-84            | Actuate to waste Beg Day<br><b>Shot weight collection</b>                                     |
|           | 85-88            | Actuate to waste Mid Day<br><b>Shot weight collection</b>                                     |
|           | 89-92            | Actuate to waste End Day<br><b>Shot weight collection</b>                                     |
| Tuesday   | 93-104           | Actuate to waste<br><b>Shot weight collection</b>   |
| Wednesday | 105-106**        | 24 hr Dose<br><b>Shot weight collection</b>   |
|           | 107-108          | Actuate to waste after 24 hr Dose<br>(actuactions 105 & 106)<br><b>Shot weight collection</b> |
|           | 109-112          | Actuate to waste Mid Day<br><b>Shot weight collection</b>                                     |
|           | 113-116          | Actuate to waste End Day<br><b>Shot weight collection</b>                                     |
| Thursday  | 117-128          | Actuate to waste<br><b>Shot weight collection</b>   |
| Friday    | 129 & 130**      | 24 hr SAR<br><b>Shot weight collection</b>  |
|           | 131-132          | Actuate to waste after 24 hr SAR<br>(actuactions 129 & 130)<br><b>Shot weight collection</b>  |
|           | 133-136          | Actuate to waste Mid Day<br><b>Shot weight collection</b>                                     |
|           | 137-140          | Actuate to waste End Day<br><b>Shot weight collection</b>                                     |

---

\*\* Must occur at least 24 hours from completion of prior wasting timepoint.

| Week 3    |                  |  |
|-----------|------------------|--|
| Day       | Actuation Number | Action   |
| Monday    | 141-144          | Actuate to waste Beg Day<br>Shot weight collection                                     |
|           | 145-148          | Actuate to waste Mid Day<br>Shot weight collection                                     |
|           | 149-152          | Actuate to waste End Day<br>Shot weight collection                                     |
| Tuesday   | 153-164          | Actuate to waste<br>Shot weight collection   |
| Wednesday | 165-166**        | 24 hr Dose<br>Shot weight collection   |
|           | 167-168          | Actuate to waste after 24 hr Dose<br>(actuactions 165 & 166)<br>Shot weight collection |
|           | 169-172          | Actuate to waste Mid Day<br>Shot weight collection                                     |
|           | 173-174          | Actuate to waste End Day<br>Shot weight collection                                     |
| Thursday  | 175-186          | Actuate to waste<br>Shot weight collection   |
| Friday    | 187 & 188**      | 24 hr SAR<br>Shot weight collection  |
|           | 189-198 *        | Clear sample inlet tube<br>(Water Content)   |
|           | 199-203 *        | Water Content Sample   |

## V. Test Methods

| Method             | ATM # | BIPI TP # | Method Transfer |
|--------------------|-------|-----------|-----------------|
| Water              | ATM-  |           | TTP-            |
| Uniformity of Dose | ATM-  |           | TTP-            |

## VI. Reporting

Report the data using the Access Database provided in a Certificate of Analysis format.  
An electronic copy of the data should also be provided.

## VII. Attachments

None.

\* May be performed as late as the Monday following the last SAR.

\*\* Must occur at least 24 hours from completion of prior wasting timepoint.

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7501

Theoretical Water Content: Inherent

Analyst: Brook White

| Actuation Number | Test | Inverted            |                    |         |       |         |       |         |       |         |       |
|------------------|------|---------------------|--------------------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1               |                    | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
|                  |      | Alb SO4             | IPB                | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 115.13              | 20.56              | 118.21  | 20.85 | 115.17  | 20.64 | 117.57  | 20.63 | 122.36  | 21.72 |
| 26               | SAR  | 113.42              | 19.74              | 112.10  | 19.45 | 117.79  | 20.54 | 117.26  | 20.34 | 119.82  | 21.42 |
| 33               | SAR  | 77.03               | 19.66              | 66.59   | 19.61 | 90.17   | 19.95 | 62.84   | 19.43 | 82.38   | 19.42 |
| 34               | SAR  | 245.57              | 20.14              | 202.70  | 19.47 | 216.22  | 19.59 | 238.14  | 21.24 | 214.94  | 20.25 |
| 45-46            | Dose | 159.74              | 19.46              | 139.27  | 18.83 | 144.63  | 19.52 | 169.23  | 19.22 | 149.47  | 19.76 |
| 69               | SAR  | 60.95               | 19.29              | 75.53   | 19.87 | 76.67   | 19.34 | 70.93   | 18.94 | 128.19  | 18.60 |
| 70               | SAR  | 240.53              | 20.28              | 237.87  | 19.34 | 256.88  | 20.88 | 232.40  | 20.84 | 222.82  | 21.00 |
| 105-106          | Dose | 133.76              | 19.27              | 95.46   | 19.54 | 95.94   | 20.03 | 92.92   | 19.78 | 95.99   | 19.87 |
| 129              | SAR  | 169.42 <sup>1</sup> | 19.62 <sup>1</sup> | 57.22   | 19.31 | 79.18   | 19.93 | 69.29   | 18.84 | 56.62   | 18.79 |
| 130              | SAR  | 161.63 <sup>2</sup> | 19.62 <sup>2</sup> | 181.90  | 19.69 | 156.13  | 20.43 | 177.50  | 20.78 | 161.25  | 19.90 |
| 165-166          | Dose | 103.76              | 19.55              | 101.53  | 19.61 | 104.51  | 20.45 | 100.96  | 19.96 | 92.65   | 19.65 |
| 187              | SAR  | 87.58               | 19.52              | 58.20   | 25.05 | 47.18   | 19.52 | 57.62   | 19.52 | 92.57   | 19.22 |
| 188              | SAR  | 146.99              | 20.16              | 168.57  | 20.82 | 142.08  | 21.13 | 153.97  | 20.93 | 136.70  | 20.66 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 263   | 257   | 260   | 262   | 265   |
| 199-203   | 364   | 343   | 343   | 331   | 325   |

<sup>1</sup>These values were confirmed by re-injection since they were not in trend with the data set. The re-injected samples yielded results of 167.16 and 18.81 µg for Albuterol Sulfate and Ipratropium Bromide respectively.

<sup>2</sup>These values were confirmed by re-injection since they were not in trend with the data set. The re-injected samples yielded results of 160.32 and 19.42 µg for Albuterol Sulfate and Ipratropium Bromide respectively.

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7501

Theoretical Water Content: Inherent

Analyst: Nichole Bynum

| Actuation Number | Test | Upright |       |                     |                    |         |       |                    |                    |                    |                    |
|------------------|------|---------|-------|---------------------|--------------------|---------|-------|--------------------|--------------------|--------------------|--------------------|
|                  |      | Can 6   |       | Can 7               |                    | Can 8   |       | Can 9              |                    | Can 10             |                    |
| 25               | SAR  | Alb SO4 | IPB   | Alb SO4             | IPB                | Alb SO4 | IPB   | Alb SO4            | IPB                | Alb SO4            | IPB                |
| 26               | SAR  | 114.77  | 20.55 | 121.29              | 21.97              | 121.85  | 21.73 | 115.84             | 20.69              | 121.53             | 21.52              |
| 33               | SAR  | 119.57  | 20.87 | 118.26              | 20.88              | 117.92  | 22.26 | 121.20             | 21.46              | 119.29             | 21.36              |
| 34               | SAR  | 39.60   | 19.58 | 47.68               | 20.48              | 42.23   | 20.70 | 40.38              | 20.39              | 47.37              | 20.29              |
| 45-46            | Dose | 164.26  | 19.83 | 155.10              | 19.75              | 156.82  | 19.18 | 169.58             | 19.66              | 156.61             | 19.84              |
| 69               | SAR  | 98.77   | 19.56 | 99.56               | 19.85              | 99.49   | 20.45 | 104.71             | 20.18              | 98.57              | 20.05              |
| 70               | SAR  | 46.50   | 19.64 | 50.50               | 19.16              | 41.45   | 19.82 | 50.31              | 20.04              | 52.55              | 20.33              |
| 105-106          | Dose | 155.93  | 19.89 | 157.65              | 20.68              | 159.42  | 19.32 | 150.98             | 18.73              | 158.94             | 20.84              |
| 129              | SAR  | 88.03   | 18.82 | 96.15               | 20.00              | 94.30   | 19.59 | 92.94              | 19.42              | 106.61             | 20.44              |
| 130              | SAR  | 51.96   | 20.29 | 55.25               | 20.24              | 77.68   | 20.06 | 89.13              | 61.90              | 78.36              | 17.70              |
| 165-166          | Dose | 150.76  | 20.17 | 141.95              | 19.57              | 142.01  | 19.74 | 130.40             | 18.61              | 136.65             | 20.80              |
| 187              | SAR  | 94.23   | 19.57 | 102.94 <sup>2</sup> | 17.66 <sup>2</sup> | 101.73  | 18.17 | 99.60 <sup>2</sup> | 16.83 <sup>2</sup> | 98.36 <sup>2</sup> | 16.97 <sup>2</sup> |
| 188              | SAR  | 48.15   | 20.44 | 60.41               | 20.65              | 42.98   | 20.22 | 46.04              | 20.68              | 60.89              | 22.05              |
|                  |      | 149.23  | 20.66 | 140.19              | 20.92              | 150.85  | 21.02 | 152.92             | 20.66              | 138.40             | 21.78              |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 6 | Can 7 | Can 8 | Can 9 | Can 10           |
|-----------|-------|-------|-------|-------|------------------|
| 16-20     | 275   | 270   | 260   | 268   | 279              |
| 199-203   | 341   | 337   | 327   | 330   | 402 <sup>1</sup> |

<sup>1</sup>The analyst observed an unusually weak plume for all 5 actuations.

<sup>2</sup>The analyst observed that the first actuation sounded weaker than normal.

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7505

Theoretical Water Content: 1500ppm

Analyst: Brook White

| Actuation Number | Test | Inverted |       |         |       |         |       |         |       |         |       |
|------------------|------|----------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 1    |       | Can 2   |       | Can 3   |       | Can 4   |       | Can 5   |       |
|                  |      | Alb SO4  | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 116.21   | 20.67 | 115.98  | 20.74 | 119.88  | 21.38 | 126.79  | 22.26 | 118.66  | 21.10 |
| 26               | SAR  | 118.46   | 20.67 | 116.09  | 20.10 | 121.20  | 21.20 | 122.00  | 21.34 | 116.72  | 20.73 |
| 33               | SAR  | 112.55   | 19.84 | 116.28  | 20.82 | 114.67  | 20.34 | 113.06  | 20.31 | 114.61  | 20.86 |
| 34               | SAR  | 117.06   | 20.50 | 111.72  | 19.32 | 118.41  | 20.54 | 110.90  | 19.19 | 115.39  | 20.39 |
| 45-46            | Dose | 110.45   | 19.80 | 111.29  | 19.78 | 114.27  | 20.55 | 112.32  | 19.91 | 113.49  | 20.17 |
| 69               | SAR  | 99.24    | 19.26 | 108.44  | 20.24 | 99.74   | 20.02 | 105.97  | 20.25 | 101.45  | 20.17 |
| 70               | SAR  | 122.24   | 20.06 | 119.12  | 20.08 | 124.01  | 20.35 | 118.92  | 19.94 | 124.88  | 20.31 |
| 105-106          | Dose | 111.73   | 20.21 | 113.12  | 20.06 | 113.20  | 20.48 | 112.03  | 19.96 | 114.97  | 20.84 |
| 129              | SAR  | 113.98   | 20.14 | 113.93  | 19.94 | 112.79  | 20.78 | 110.89  | 20.07 | 108.77  | 19.69 |
| 130              | SAR  | 118.36   | 21.38 | 115.40  | 20.20 | 116.83  | 20.81 | 115.00  | 20.45 | 125.57  | 22.64 |
| 165-166          | Dose | 107.47   | 19.95 | 112.14  | 20.24 | 110.44  | 20.25 | 109.26  | 20.24 | 109.89  | 20.29 |
| 187              | SAR  | 117.19   | 20.72 | 119.51  | 20.61 | 110.54  | 20.51 | 107.43  | 20.37 | 119.79  | 20.74 |
| 188              | SAR  | 107.92   | 20.35 | 114.33  | 20.87 | 122.09  | 22.64 | 112.39  | 21.22 | 117.01  | 22.25 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 1 | Can 2 | Can 3 | Can 4 | Can 5 |
|-----------|-------|-------|-------|-------|-------|
| 16-20     | 1381  | 1430  | 1361  | 1413  | 1380  |
| 199-203   | 1488  | 1488  | 1438  | 1503  | 1486  |

# Uniformity of Dose Water Challenge

Results reported in µg recovered

Lot Number:

AS Number: 02-10-7505

Theoretical Water Content: 1500ppm

Analyst: Nichole Bynum

| Actuation Number | Test | Upright |       |         |       |         |       |         |       |         |       |
|------------------|------|---------|-------|---------|-------|---------|-------|---------|-------|---------|-------|
|                  |      | Can 6   |       | Can 7   |       | Can 8   |       | Can 9   |       | Can 10  |       |
|                  |      | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   | Alb SO4 | IPB   |
| 25               | SAR  | 116.73  | 20.84 | 117.75  | 20.66 | 119.88  | 20.65 | 117.66  | 20.83 | 125.01  | 21.96 |
| 26               | SAR  | 113.69  | 20.06 | 115.37  | 20.18 | 118.13  | 20.57 | 125.89  | 21.77 | 124.74  | 21.81 |
| 33               | SAR  | 119.64  | 22.47 | 116.18  | 19.72 | 118.66  | 20.09 | 119.61  | 20.44 | 116.61  | 20.40 |
| 34               | SAR  | 119.72  | 19.94 | 117.46  | 19.94 | 118.13  | 19.56 | 119.28  | 19.61 | 121.44  | 20.13 |
| 45-46            | Dose | 114.78  | 19.30 | 114.70  | 19.16 | 108.98  | 18.65 | 114.77  | 19.74 | 112.26  | 18.47 |
| 69               | SAR  | 112.18  | 20.14 | 116.05  | 19.46 | 111.63  | 19.43 | 114.12  | 20.20 | 114.19  | 19.60 |
| 70               | SAR  | 122.24  | 20.85 | 117.55  | 19.33 | 110.99  | 18.25 | 118.92  | 20.47 | 118.45  | 18.81 |
| 105-106          | Dose | 116.16  | 19.95 | 114.83  | 19.56 | 112.33  | 19.05 | 112.30  | 20.42 | 114.25  | 19.85 |
| 129              | SAR  | 113.66  | 20.70 | 109.67  | 19.65 | 122.59  | 19.86 | 116.34  | 20.45 | 112.39  | 19.84 |
| 130              | SAR  | 113.45  | 20.18 | 121.31  | 19.92 | 109.40  | 18.51 | 112.48  | 19.30 | 121.42  | 20.02 |
| 165-166          | Dose | 109.10  | 20.13 | 108.94  | 19.64 | 114.02  | 18.91 | 110.36  | 20.11 | 115.55  | 17.60 |
| 187              | SAR  | 97.28   | 20.74 | 99.90   | 19.93 | 105.47  | 20.31 | 107.48  | 20.96 | 104.49  | 20.53 |
| 188              | SAR  | 126.94  | 21.01 | 116.29  | 19.90 | 111.20  | 19.25 | 111.98  | 20.10 | 110.54  | 20.38 |

Alb SO4 = Albuterol Sulfate

IPB = Ipratropium Bromide

SAR = Single Actuation Reproducibility

## Water Content (ppm)

| Actuation | Can 6 | Can 7 | Can 8 | Can 9 | Can 10 |
|-----------|-------|-------|-------|-------|--------|
| 16-20     | 1367  | 1393  | 1369  | 1428  | 1381   |
| 199-203   | 1491  | 1492  | 1471  | 1562  | 1498   |



**Boehringer  
Ingelheim**

**ANALYTICAL SCIENCES DEPARTMENT**

**ANALYTICAL REPORT**

**Title:** Results for In Use 6 Hour Delay Single Actuation Reproducibility Study for Ipratropium Bromide/Albuterol Sulfate Inhalation Aerosol, with Water Contents Ranging from Inherent Water to 2500 ppm

**Author(s):** Johanna Ubben

**Date of Report:** 05 Feb 2003

**Product / Project:**

**Project Code:** BP

**Requester:** George DeStefano

**cc:** P. D. Curry, Jr, Ph.D.; E. Gump, Ph.D.; R. Patel, Ph.D.

**Sample ID:** 02-10-7515, 02-10-7516, 02-10-7517, 02-10-7518, 02-10,7519, 02-10-7520

**Key Words:** Combivent HFA, dose uniformity, water content, single actuation reproducibility, valve delivery through life, total can assay

**File Number:** AR-030012

**Security:** Read Only

**Signature:** *Johanna Ubben* Feb 21, 2003

**Verification:** *[Signature]* Feb. 21, 2003

**Data Reference(s):** 5886/083

## **I. Introduction**

Samples from Ipratropium Bromide/Albuterol Sulfate Inhalation Aerosol, Lot: \_\_\_\_\_ were submitted to the Pulmonary Analysis Laboratory for a patient in-use, single actuation reproducibility study. The protocol detailing the experimental plan for this study is shown in Attachment 1. The purpose of this report is to summarize the results of this testing

## **II. Samples**

All samples used the Bepak Valve #BK0080332, 0.05% PVP, and 10 % ethanol.

- \_\_\_\_\_ A1 – Inherent Water
- \_\_\_\_\_ A2 – 800 ppm Water
- \_\_\_\_\_ A3 – 1000 ppm Water
- \_\_\_\_\_ A4 – 1200 ppm Water
- \_\_\_\_\_ A5 – 1500 ppm Water
- \_\_\_\_\_ A6 – 2500 ppm Water

## **III. Test Procedures**

TP-\_\_\_\_\_ Water Content

TP-\_\_\_\_\_ Dose Uniformity (US)

TP-\_\_\_\_\_ Total Can Assay

Single Actuation Reproducibility – Experimental Procedure described in protocol

TP-\_\_\_\_\_, Valve Delivery Through Life – with modification described in protocol



|         | Orientation | Canister<br>Number | Beginning<br>of Study<br>(ppm) | End of<br>Study<br>(ppm) |
|---------|-------------|--------------------|--------------------------------|--------------------------|
| 10/1/94 | Inverted    | 3                  | 1388.0                         | 1459.0                   |
| 10/1/94 | Inverted    | 4                  | 1383.0                         | 1612.0                   |
| 10/1/94 | Inverted    | 5                  | 1375.0                         | 1597.0                   |
| 10/1/94 | Upright     | 1                  | 2724.0                         | 1799.0                   |
| 10/1/94 | Upright     | 2                  | 4774.0                         | 1767.0                   |
| 10/1/94 | Upright     | 3                  | 2347.0                         | 1791.0                   |
| 10/1/94 | Upright     | 4                  | 2318.0                         | 1840.0                   |
| 10/1/94 | Upright     | 5                  | 2277.0                         | 1847.0                   |
| 10/1/94 | Inverted    | 1                  | 2276.0                         | 1877.0                   |
| 10/1/94 | Inverted    | 2                  | 2365.0                         | 2571.0                   |
| 10/1/94 | Inverted    | 3                  | 2313.0                         | 1781.0                   |
| 10/1/94 | Inverted    | 4                  | 2333.0                         | 2002.0                   |
| 10/1/94 | Inverted    | 5                  | 2344.0                         | 2572.0                   |

#### B. Single Actuation Reproducibility after 6 Hour Delay

Ten canisters from each lot were tested for single actuation reproducibility following a 6 hour rest period after priming the canisters. Five canisters from each lot were stored upright and five canisters were stored inverted. The results for each lot are shown in Tables 2, 3, 4, 5, 6, and 7. The albuterol sulfate data from these tables was also compiled and is presented graphically in Figure 2.

For lots 10/1/94, 10/2/94, and 10/3/94 testing was discontinued after the first week. These lots exhibited the albuterol sulfate trend of being low for the first actuation 1 and high for the second actuation with a 6 hour rest after priming. The remaining lots do not exhibit this behavior and were tested through the life of the canister following the patient in use actuation protocol.

Table 2. Single Actuation Reproducibility for Lot Upright and Inverted Samples, Inherent Water

| Albuterol Sulfate    |                                 |                                 |                                 |                                 |                                 |                                 |
|----------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| Lot /Orientation/Can | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
|                      | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 1 -Up-1          | 49.87                           | 126.95                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Up-2          | 53.22                           | 127.20                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Up-3          | 56.08                           | 124.65                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Up-4          | 51.62                           | 127.27                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Up-5          | 50.58                           | 127.07                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-1         | 78.10                           | 151.04                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-2         | 113.18                          | 154.54                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-3         | 88.07                           | 174.11                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-4         | 93.35                           | 169.84                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-5         | 70.53                           | 161.62                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Ipratropium Bromide  |                                 |                                 |                                 |                                 |                                 |                                 |
| Lot /Orientation/Can | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
|                      | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 2 -Up-1          | 95.24                           | 98.92                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-2          | 93.13                           | 104.19                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-3          | 99.05                           | 102.11                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-4          | 96.32                           | 101.88                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-5          | 98.14                           | 96.52                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-1         | 92.54                           | 95.59                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-2         | 94.69                           | 91.53                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-3         | 89.15                           | 100.92                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-4         | 93.30                           | 93.63                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-5         | 96.02                           | 101.16                          | N/A                             | N/A                             | N/A                             | N/A                             |

Table 3. Single Actuation Reproducibility for Lot 1  
Upright and Inverted Samples, 800 ppm Water

| Albuterol Sulfate    |                                 |                                 |                                 |                                 |                                 |                                 |
|----------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| Lot /Orientation/Can | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
|                      | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 1 -Up-1          | 48.60                           | 137.69                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Up-2          | 49.35                           | 123.03                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Up-3          | 59.65                           | 123.56                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Up-4          | 49.85                           | 128.45                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Up-5          | 45.63                           | 128.28                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-1         | 89.50                           | 137.83                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-2         | 63.97                           | 145.05                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-3         | 71.63                           | 122.67                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-4         | 75.94                           | 126.01                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 1 -Inv-5         | 65.85                           | 137.82                          | N/A                             | N/A                             | N/A                             | N/A                             |

| Ipratropium Bromide  |                                 |                                 |                                 |                                 |                                 |                                 |
|----------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| Lot /Orientation/Can | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
|                      | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 2 -Up-1          | 98.15                           | 99.55                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-2          | 97.86                           | 96.72                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-3          | 95.91                           | 103.24                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-4          | 98.77                           | 97.53                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-5          | 96.66                           | 101.29                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-1         | 95.75                           | 98.33                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-2         | 93.67                           | 99.81                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-3         | 97.13                           | 97.87                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-4         | 95.95                           | 97.28                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-5         | 95.27                           | 102.27                          | N/A                             | N/A                             | N/A                             | N/A                             |

Table 4. Single Actuation Reproducibility for Lot 2  
Upright and Inverted Samples, 1000 ppm Water

| Albuterol Sulfate    |                                 |                                 |                                 |                                 |                                 |                                 |
|----------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| Lot /Orientation/Can | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
|                      | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 2 -Up-1          | 71.63                           | 110.75                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-2          | 56.02                           | 130.26                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-3          | 51.97                           | 126.82                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-4          | 52.41                           | 124.41                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-5          | 80.63                           | 104.43                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-1         | 63.38                           | 128.71                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-2         | 65.70                           | 128.20                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-3         | 73.37                           | 118.64                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-4         | 77.65                           | 121.92                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-5         | 66.54                           | 126.07                          | N/A                             | N/A                             | N/A                             | N/A                             |

| Ipratropium Bromide  |                                 |                                 |                                 |                                 |                                 |                                 |
|----------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| Lot /Orientation/Can | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
|                      | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 2 -Up-1          | 95.35                           | 96.41                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-2          | 96.43                           | 99.95                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-3          | 97.69                           | 97.97                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-4          | 93.75                           | 99.16                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Up-5          | 95.13                           | 96.11                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-1         | 95.56                           | 100.24                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-2         | 96.95                           | 100.28                          | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-3         | 97.81                           | 97.64                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-4         | 97.71                           | 98.70                           | N/A                             | N/A                             | N/A                             | N/A                             |
| Lot 2 -Inv-5         | 96.17                           | 100.76                          | N/A                             | N/A                             | N/A                             | N/A                             |

**Table 5. Single Actuation Reproducibility for Lot Upright and Inverted Samples, 1200 ppm Water**

| Albuterol Sulfate  |  |             |             |             |             |             |             |
|--|--|-------------|-------------|-------------|-------------|-------------|-------------|
|  |  | Week 1      |             | Week 2      |             | Week 3      |             |
|  |  | Actuation 1 | Actuation 2 | Actuation 1 | Actuation 2 | Actuation 1 | Actuation 2 |
|  |  | (% of       | (% of       | (% of       | (% of       | (% of       | (% of       |
| Lot /Orientation/Can   |  | Theory)     | Theory)     | Theory)     | Theory)     | Theory)     | Theory)     |
| 00 |  |             |             |             |             |             |             |

Table 6. Single Actuation Reproducibility for Lot 1  
Upright and Inverted Samples, 1500 ppm Water

| Albuterol Sulfate    |                                 |                                 |                                 |                                 |                                 |                                 |
|----------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| Lot /Orientation/Can | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
|                      | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 1 -Up-1          | 94.62                           | 99.02                           | 102.73                          | 110.10                          | 101.00                          | 108.91                          |
| Lot 1 -Up-2          | 92.30                           | 101.31                          | 100.62                          | 105.94                          | 91.17                           | 112.39                          |
| Lot 1 -Up-3          | 94.41                           | 101.03                          | 99.17                           | 116.08                          | 101.08                          | 107.24                          |
| Lot 1 -Up-4          | 95.39                           | 99.19                           | 95.33                           | 105.99                          | 92.86                           | 109.29                          |
| Lot 1 -Up-5          | 98.04                           | 100.37                          | 95.85                           | 110.52                          | 89.44                           | 110.49                          |
| Lot 1 -Inv-1         | 96.36                           | 99.38                           | 92.49                           | 109.67                          | 83.09                           | 107.79                          |
| Lot 1 -Inv-2         | 94.78                           | 105.63                          | 93.82                           | 116.65                          | 85.19                           | 117.40                          |
| Lot 1 -Inv-3         | 98.08                           | 101.14                          | 101.80                          | 107.13                          | 94.68                           | 110.60                          |
| Lot 1 -Inv-4         | 94.92                           | 101.36                          | 94.04                           | 108.36                          | 88.57                           | 108.73                          |
| Lot 1 -Inv-5         | 93.20                           | 99.04                           | 99.38                           | 110.70                          | 94.67                           | 102.19                          |

| Ipratropium Bromide  |                                 |                                 |                                 |                                 |                                 |                                 |
|----------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| Lot /Orientation/Can | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
|                      | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 2 -Up-1          | 97.57                           | 97.84                           | 95.93                           | 100.56                          | 104.01                          | 105.51                          |
| Lot 2 -Up-2          | 96.13                           | 99.37                           | 99.04                           | 97.10                           | 102.07                          | 107.00                          |
| Lot 2 -Up-3          | 101.21                          | 98.24                           | 99.87                           | 102.88                          | 107.59                          | 107.83                          |
| Lot 2 -Up-4          | 98.40                           | 98.62                           | 98.13                           | 97.41                           | 101.62                          | 104.50                          |
| Lot 2 -Up-5          | 100.59                          | 98.85                           | 97.67                           | 101.23                          | 101.67                          | 106.88                          |
| Lot 2 -Inv-1         | 97.29                           | 98.37                           | 94.27                           | 102.27                          | 91.78                           | 105.21                          |
| Lot 2 -Inv-2         | 97.21                           | 102.94                          | 98.97                           | 105.40                          | 103.99                          | 108.79                          |
| Lot 2 -Inv-3         | 100.14                          | 100.35                          | 97.97                           | 101.26                          | 102.20                          | 107.54                          |
| Lot 2 -Inv-4         | 97.44                           | 101.24                          | 98.73                           | 100.17                          | 103.00                          | 107.15                          |
| Lot 2 -Inv-5         | 95.39                           | 97.97                           | 95.20                           | 101.44                          | 99.90                           | 106.96                          |

Table 7. Single Actuation Reproducibility for Lot Upright and Inverted Samples, 2500 ppm Water

| Albuterol Sulfate    |        |                                 |                                 |                                 |                                 |                                 |                                 |
|----------------------|--------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
|                      |        | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
| Lot /Orientation/Can |        | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 1                | -Up-1  | 95.51                           | 96.18                           | 93.06                           | 106.82                          | 89.02                           | 98.71                           |
|                      | -Up-2  | 91.40                           | 98.26                           | 97.25                           | 100.72                          | 91.40                           | 96.23                           |
|                      | -Up-3  | 97.37                           | 98.65                           | 99.53                           | 96.12                           | 90.59                           | 103.09                          |
|                      | -Up-4  | 95.08                           | 95.98                           | 91.04                           | 103.12                          | 79.05                           | 103.87                          |
|                      | -Up-5  | 95.89                           | 100.90                          | 104.42                          | 95.22                           | 92.35                           | 100.06                          |
| Lot 2                | -Inv-1 | 99.18                           | 97.61                           | 98.05                           | 100.63                          | 89.50                           | 104.01                          |
|                      | -Inv-2 | 98.00                           | 99.87                           | 103.83                          | 105.87                          | 87.10                           | 97.13                           |
|                      | -Inv-3 | 95.97                           | 98.53                           | 99.28                           | 103.49                          | 87.53                           | 96.46                           |
|                      | -Inv-4 | 97.42                           | 97.94                           | 95.65                           | 102.85                          | 88.18                           | 97.42                           |
|                      | -Inv-5 | 100.73                          | 103.16                          | 103.08                          | 105.17                          | 91.20                           | 100.93                          |
| Ipratropium Bromide  |        |                                 |                                 |                                 |                                 |                                 |                                 |
|                      |        | Week 1                          |                                 | Week 2                          |                                 | Week 3                          |                                 |
| Lot /Orientation/Can |        | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) | Actuation 1<br>(% of<br>Theory) | Actuation 2<br>(% of<br>Theory) |
| Lot 1                | -Up-1  | 98.49                           | 98.44                           | 99.09                           | 101.27                          | 103.15                          | 103.73                          |
|                      | -Up-2  | 92.51                           | 100.17                          | 99.39                           | 95.11                           | 103.77                          | 101.20                          |
|                      | -Up-3  | 99.17                           | 98.09                           | 99.57                           | 95.75                           | 100.00                          | 105.30                          |
|                      | -Up-4  | 97.77                           | 96.62                           | 95.22                           | 97.40                           | 91.87                           | 109.92                          |
|                      | -Up-5  | 95.55                           | 98.96                           | 99.36                           | 95.79                           | 102.33                          | 100.74                          |
| Lot 2                | -Inv-1 | 98.09                           | 99.33                           | 97.88                           | 100.19                          | 102.34                          | 110.94                          |
|                      | -Inv-2 | 98.69                           | 99.38                           | 103.30                          | 106.59                          | 102.82                          | 109.17                          |
|                      | -Inv-3 | 94.24                           | 101.02                          | 96.22                           | 103.06                          | 100.14                          | 108.23                          |
|                      | -Inv-4 | 97.41                           | 96.45                           | 96.90                           | 100.47                          | 102.63                          | 107.16                          |
|                      | -Inv-5 | 96.72                           | 102.32                          | 102.22                          | 105.63                          | 104.23                          | 110.85                          |

**C. Dose Uniformity after 6 Hour Delay**

Ten canisters from each lot were tested for uniformity of dose (2 actuations) following a 6 hour rest period after priming the canisters. Five canisters from each lot were stored upright and five canisters were stored inverted. The results for each lot are shown in Tables 8, 9, 10, 11, 12, and 13. The albuterol sulfate data from these tables was also compiled and is presented graphically in Figures 3.

For lots testing was discontinued after the first week. These lots had exhibited poor single actuation reproducibility during the first week of testing. The remaining lots were tested through the life of the canister following the patient in use actuation protocol.



Table 8. Uniformity of Dose for Lot (Inherent Water), Upright and Inverted Samples

| Lot/Orientation/Can | Week 1                             |   |  | Week 2                             |   |  | Week 3                             |   |  |
|---------------------|------------------------------------|---|--|------------------------------------|---|--|------------------------------------|---|--|
|                     | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) |  | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) |  | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) |  |
| Up/1                | 93.58                              | 99.67                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Up/2                | 91.59                              | 100.38                                  |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Up/3                | 90.11                              | 98.39                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Up/4                | 87.17                              | 95.25                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Up/5                | 95.20                              | 99.56                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/1               | 148.63                             | 97.29                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/2               | 153.88                             | 97.25                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/3               | 168.93                             | 102.22                                  |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/4               | 125.87                             | 96.76                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/5               | 118.27                             | 99.42                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |

Table 9. Uniformity of Dose for Lot: (800 ppm Water), Upright and Inverted Samples

| Lot/Orientation/Can | Week 1                             |   |  | Week 2                             |   |  | Week 3                             |   |  |
|---------------------|------------------------------------|---|--|------------------------------------|---|--|------------------------------------|---|--|
|                     | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) |  | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) |  | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) |  |
| /Up/1               | 100.61                             | 98.52                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| /Up/2               | 88.43                              | 98.38                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| /Up/3               | 95.12                              | 98.97                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| /Up/4               | 89.92                              | 96.90                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| /Up/5               | 88.60                              | 97.17                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| /Inv/1              | 101.06                             | 97.01                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| /Inv/2              | 102.17                             | 97.76                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| /Inv/3              | 101.22                             | 99.41                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| /Inv/4              | 105.29                             | 97.20                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| /Inv/5              | 99.46                              | 99.34                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |

Table 10. Uniformity of Dose for Lot (1000 ppm Water), Upright and Inverted Samples

| Lot/Orientation/Can | Week 1                             |   |  | Week 2                             |   |  | Week 3                             |   |  |
|---------------------|------------------------------------|---|--|------------------------------------|---|--|------------------------------------|---|--|
|                     | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) |  | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) |  | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) |  |
| Up/1                | 94.82                              | 101.07                                  |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Up/2                | 92.18                              | 97.01                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Up/3                | 94.90                              | 102.35                                  |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Up/4                | 87.00                              | 94.65                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Up/5                | 91.39                              | 97.60                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/1               | 103.68                             | 100.80                                  |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/2               | 101.18                             | 100.45                                  |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/3               | 101.69                             | 100.84                                  |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/4               | 97.50                              | 98.91                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |
| Inv/5               | 99.93                              | 99.46                                   |  | N/A                                | N/A                                     |  | N/A                                | N/A                                     |  |

Table 11. Uniformity of Dose for Lot (1200 ppm Water), Upright and Inverted Samples

| Lot/Orientation/Can | Week 1            |               |                                   | Week 2            |               |                                   | Week 3            |               |                                   |
|---------------------|-------------------|---------------|-----------------------------------|-------------------|---------------|-----------------------------------|-------------------|---------------|-----------------------------------|
|                     | Albuterol Sulfate |               | Ipratropium Bromide (% of Theory) | Albuterol Sulfate |               | Ipratropium Bromide (% of Theory) | Albuterol Sulfate |               | Ipratropium Bromide (% of Theory) |
|                     | (% of Theory)     | (% of Theory) |                                   | (% of Theory)     | (% of Theory) |                                   | (% of Theory)     | (% of Theory) |                                   |
| Up/1                | 98.05             | 96.57         | 96.57                             | 94.30             | 98.98         | 98.98                             | 101.45            | 101.40        |                                   |
| Up/2                | 100.06            | 98.15         | 98.15                             | 100.77            | 99.43         | 99.43                             | 103.26            | 99.80         |                                   |
| Up/3                | 94.51             | 97.32         | 97.32                             | 83.48             | 99.77         | 99.77                             | 103.59            | 101.52        |                                   |
| Up/4                | 97.56             | 97.49         | 97.49                             | 96.82             | 99.70         | 99.70                             | 100.46            | 100.63        |                                   |
| Up/5                | 89.82             | 99.02         | 99.02                             | 93.53             | 104.60        | 104.60                            | 108.98            | 106.29        |                                   |
| Inv/1               | 99.90             | 100.52        | 100.52                            | 97.95             | 101.21        | 101.21                            | 103.23            | 104.73        |                                   |
| Inv/2               | 103.60            | 102.42        | 102.42                            | 102.00            | 104.31        | 104.31                            | 100.12            | 101.44        |                                   |
| Inv/3               | 94.77             | 92.73         | 92.73                             | 97.57             | 100.38        | 100.38                            | 97.78             | 102.36        |                                   |
| Inv/4               | 99.66             | 98.75         | 98.75                             | 96.92             | 102.10        | 102.10                            | 99.31             | 103.07        |                                   |
| Inv/5               | 102.07            | 99.49         | 99.49                             | 95.21             | 100.01        | 100.01                            | 96.91             | 104.22        |                                   |

Table 12. Uniformity of Dose for Lot (1500 ppm Water), Upright and Inverted Samples

| Lot/Orientation/Can | Week 1                             |   |                                    | Week 2                                  |                                    |   | Week 3                             |   |                                    |
|---------------------|------------------------------------|---|------------------------------------|---|------------------------------------|---|------------------------------------|---|------------------------------------|
|                     | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) | Albuterol Sulfate<br>(% of Theory) | Ipratropium<br>Bromide (% of<br>Theory) | Albuterol Sulfate<br>(% of Theory) |
| Upr/1               | 101.37                             | 100.07                                  | 101.25                             | 100.79                                  | 103.62                             | 102.01                                  |                                    |   |                                    |
| Upr/2               | 99.40                              | 97.58                                   | 102.79                             | 100.28                                  | 99.08                              | 101.35                                  |                                    |   |                                    |
| Upr/3               | 102.85                             | 100.26                                  | 103.15                             | 101.15                                  | 99.52                              | 101.78                                  |                                    |   |                                    |
| Upr/4               | 102.06                             | 99.23                                   | 102.39                             | 101.28                                  | 98.10                              | 97.21                                   |                                    |   |                                    |
| Upr/5               | 102.73                             | 100.78                                  | 102.55                             | 100.73                                  | 104.53                             | 102.70                                  |                                    |   |                                    |
| Inv/1               | 99.48                              | 99.40                                   | 98.77                              | 100.67                                  | 100.56                             | 101.62                                  |                                    |   |                                    |
| Inv/2               | 101.03                             | 101.57                                  | 100.07                             | 102.87                                  | 98.95                              | 103.23                                  |                                    |   |                                    |
| Inv/3               | 103.35                             | 101.23                                  | 100.68                             | 101.15                                  | 102.45                             | 101.75                                  |                                    |   |                                    |
| Inv/4               | 100.28                             | 101.19                                  | 101.73                             | 102.49                                  | 99.89                              | 103.29                                  |                                    |   |                                    |
| Inv/5               | 99.35                              | 99.52                                   | 100.22                             | 101.33                                  | 99.13                              | 102.28                                  |                                    |   |                                    |

Table 13. Uniformity of Dose for Lot (2500 ppm Water), Upright and Inverted Samples

| Lot/Orientation/Can | Week 1            |               |                                   | Week 2            |               |                                   | Week 3            |               |                                   |
|---------------------|-------------------|---------------|-----------------------------------|-------------------|---------------|-----------------------------------|-------------------|---------------|-----------------------------------|
|                     | Albuterol Sulfate |               | Ipratropium Bromide (% of Theory) | Albuterol Sulfate |               | Ipratropium Bromide (% of Theory) | Albuterol Sulfate |               | Ipratropium Bromide (% of Theory) |
|                     | (% of Theory)     | (% of Theory) |                                   | (% of Theory)     | (% of Theory) |                                   | (% of Theory)     | (% of Theory) |                                   |
| Up/1                | 96.86             | 97.73         | 97.73                             | 97.96             | 97.96         | 100.81                            | 94.23             | 97.83         | 97.83                             |
| Up/2                | 102.17            | 98.49         | 98.49                             | 95.56             | 95.56         | 98.81                             | 91.56             | 91.38         | 91.38                             |
| Up/3                | 100.03            | 98.82         | 98.82                             | 97.63             | 97.63         | 97.85                             | 107.24            | 104.78        | 104.78                            |
| Up/4                | 100.72            | 102.19        | 102.19                            | 97.81             | 97.81         | 100.22                            | 89.04             | 92.07         | 92.07                             |
| Up/5                | 100.05            | 98.69         | 98.69                             | 100.00            | 100.00        | 99.72                             | 97.59             | 93.49         | 93.49                             |
| Inv/1               | 99.99             | 99.29         | 99.29                             | 96.75             | 96.75         | 101.20                            | 92.36             | 99.61         | 99.61                             |
| Inv/2               | 100.29            | 99.93         | 99.93                             | 95.69             | 95.69         | 101.00                            | 91.93             | 101.48        | 101.48                            |
| Inv/3               | 100.94            | 99.01         | 99.01                             | 97.78             | 97.78         | 100.57                            | 92.58             | 100.08        | 100.08                            |
| Inv/4               | 100.45            | 98.31         | 98.31                             | 97.59             | 97.59         | 100.65                            | 94.01             | 101.29        | 101.29                            |
| Inv/5               | 104.03            | 101.43        | 101.43                            | 101.26            | 101.26        | 103.90                            | 96.34             | 103.60        | 103.60                            |

## V. Conclusions

Water contents for all samples were in the expected ranges, although the measured results tended to be 100 – 150 ppm below target.

Single actuation results for lots with approximately 300 to 1000 ppm exhibited poor reproducibility between actuations at the beginning of the study. Testing on these canisters was discontinued after the first week of the study. Samples with 1200 ppm of water had only one can at the beginning of the study with poor reproducibility. Additional canisters from this lot started to exhibit poor reproducibility during the third week of the study. Samples with 1500 ppm water and 2500 ppm water exhibited good reproducibility between actuations throughout the study.

Dose uniformity results for all samples were within the expected ranges, with the exception of the inherent water inverted samples. These samples had 4 canisters out of 5 in excess of 125% of theory for drug delivery through the valve after a 6 hour pause. This trend was also observed for these canisters during single actuation reproducibility testing as well.

Total can assay results for the lots tested were within the expected range, with the exception of one inherent water inverted sample, which sprayed during sample preparation, resulting in low results for both albuterol sulfate and ipratropium bromide. This testing was added after the protocol was executed to confirm the analyte concentrations of the bulk solution in the canisters and that there were no losses in the metering chamber during actuation. There is good agreement of these results with the measured dose delivered through the valve, confirming that there is no loss occurring in the metering chamber during actuation.

---

The valve delivery results for all canisters tested were within expected ranges, with one exception for a canister near the end of the study, which was high and could not be explained. Since the value was high, it could not have been the result of a sticking valve, which produces low valve deliveries.






Figure 2: Single Actuation Reproducibility with 6 Hour Delay In Use – Albuterol Sulfate Results




Figure 3: Uniformity of Dose with 6 Hour Delay In Use Study – Albuterol Sulfate Results

**Attachment 1**

## Bespak 6 Hour Single Actuation and Dose through Life

### Purpose:

The purpose of this protocol is to study the effect of water on single actuation reproducibility (SAR) and Dose through life (ca. 200 actuations). Samples were specially made to use the same concentrate, but have different water contents.

### Samples:

| AS Testing Number | Lot Number | Target ppm Water |
|-------------------|------------|------------------|
| AS 02-10-7515     |            | Inherent         |
| AS 02-10-7516     |            | 800              |
| AS 02-10-7517     |            | 1000             |
| AS 02-10-7518     |            | 1200             |
| AS 02-10-7519     |            | 1500             |
| AS 02-10-7520     |            | 2500             |

Five samples in each orientation (upright & Inverted) will be tested. This will result in a total of 60 canisters to be tested. Each canister will be tested according to the sampling protocol provided.

For all methods, calculate sample weights as follows:  **$W2-W1$  = Sample Weight**

W2 is the canister weight recorded immediately following each test point

W1 is the canister weight prior to testing (previous time point's W2.)

Wasting actuations will be performed using the alternate mouthpiece design. A delay of no less than 2 hours must occur between completion of previous wasting time point and

initiation of next wasting time point when multiple wasting sessions are required on the same day. Firing and shaking to be performed as outlined in the Uniformity of Dose procedure for wasting canisters. Do not clean the interior of the valve stem following wasting. Wasting will be performed per the provided sampling protocol.

Samples must equilibrate for 6 hours before any Uniformity of Dose sampling can be performed. The single actuation reproducibility protocol (SAR) requires that two single actuations be collected in separate containers where indicated. Do not shake canisters between the individual actuations for the SAR protocol. Uniformity of Dose samples will be collected per the provided sampling protocol.

If canisters containing lower levels of water (<1100 ppm) fail SAR during the first week of testing, then they will not be tested for the remainder of the study. Zero hour testing may be used to ensure that SAR failures are truly performance related.

#### Test Methods:

(Water) As per the sampling protocol below.

(Valve Delivery) With the exception that multiple actuations will be used for the determinations and multiple samples will be collected for each canister according to the sampling plan below. Weights will be recorded at the beginning of the study and after each actuation sequence

(Dose Uniformity) With the exception that samples will be collected after allowing the canisters to sit for 6 hours. Also samples will be collected using a single actuation reproducibility (SAR) protocol, which calls for allowing the sample to sit for 6 hours then collecting two single actuations in separate containers. Samples will be taken per the protocol below.

| Start     |                  |                             |
|-----------|------------------|-----------------------------|
| Day       | Actuation Number | Action                      |
| Friday    | 1-5              | Priming                     |
|           | 6-15             | Clearing the Inlet tube     |
|           | 16-20            | Water Sample                |
| Monday    | 21-24            | Actuate to Waste Beg Day    |
|           | 25&26            | 6 hr SAR *                  |
|           | 27-32            | Actuate to Waste after SAR  |
| Tuesday   | 33-36            | Actuate to Waste Beg Day    |
|           | 37-40            | Actuate to Waste Mid Day    |
|           | 41-44            | Actuate to Waste End Day    |
| Wednesday | 45-48            | Actuate to Waste Beg Day    |
|           | 49&50            | 6 hr Dose *                 |
|           | 51-56            | Actuate to Waste after Dose |
| Thursday  | 57-60            | Actuate to Waste Beg Day    |
|           | 61-64            | Actuate to Waste Mid Day    |
|           | 65-68            | Actuate to Waste End Day    |
| Friday    | 69-72            | Actuate to Waste Beg Day    |
|           | 73-76            | Actuate to Waste Mid Day    |
|           | 77-80            | Actuate to Waste End Day    |

\* 6 hrs after actuating to waste

| Middle    |                  |                             |
|-----------|------------------|-----------------------------|
| Day       | Actuation Number | Action                      |
| Monday    | 81-84            | Actuate to Waste Beg Day    |
|           | 85-88            | Actuate to Waste Mid Day    |
|           | 89-92            | Actuate to Waste End Day    |
| Tuesday   | 93-96            | Actuate to Waste Beg Day    |
|           | 97&98            | 6 hr SAR *                  |
|           | 99-104           | Actuate to Waste after SAR  |
| Wednesday | 105-108          | Actuate to Waste Beg Day    |
|           | 109&112          | Actuate to Waste Mid Day    |
|           | 113-116          | Actuate to Waste End Day    |
| Thursday  | 117-120          | Actuate to Waste Beg Day    |
|           | 121&122          | 6 hr Dose *                 |
|           | 123-128          | Actuate to Waste after Dose |
| Friday    | 129-132          | Actuate to Waste Beg Day    |
|           | 133-136          | Actuate to Waste Mid Day    |
|           | 137-140          | Actuate to Waste after SAR  |

\* 6 hrs after actuating to waste

| End       |                  |                             |
|-----------|------------------|-----------------------------|
| Day       | Actuation Number | Action                      |
| Monday    | 141-144          | Actuate to Waste Beg Day    |
|           | 145-148          | Actuate to Waste Mid Day    |
|           | 149-152          | Actuate to Waste End Day    |
| Tuesday   | 153-156          | Actuate to Waste Beg Day    |
|           | 157-160          | Actuate to Waste Mid Day    |
|           | 161-164          | Actuate to Waste End Day    |
| Wednesday | 165-168          | Actuate to Waste Beg Day    |
|           | 169&170          | 6 hr Dose *                 |
|           | 171-174          | Actuate to Waste after Dose |
| Thursday  | 175-178          | Actuate to Waste Beg Day    |
|           | 179-182          | Actuate to Waste Mid Day    |
|           | 183-186          | Actuate to Waste End Day    |
| Friday    | 187-190          | Actuate to Waste Beg Day    |
|           | 191&192          | 6 hr SAR *                  |
| Monday    | 193-202          | Clear sample inlet tube     |
|           | 203-207          | Water – TP-00479            |

\* 6 hrs after actuating to waste

Data Reporting:

Report the data using the Access Database provided.